

harvesthealthcare[®]

General User/ Safety Guide ROYAL/ ROYAL EXTRA REPLACEMENT SYSTEM



ROYA

••••

www.harvesthealthcare.co.uk

CONTENTS

WAR	NINGS & CAUTIONS	4-
GENE	ERAL INFORMATION	6-10
1	DEFINITION OF THE GROUPS MENTIONED	7
2	NON-COMPLIANT USE	7
3	SAFETY INSTRUCTIONS	7
3.1	GENERAL SAFETY INSTRUCTIONS	7
3.2	SAFETY INFORMATION FOR THE OPERATOR	8
3.3	SAFETY INFORMATION FOR THE USER	8
3.4	SYMBOLS USED	9
3.5	CLEANING & DISINFECTION	10
3.6	SERVICING & MAINTENANCE	10
3.7	SERVICE LIFE & DISPOSAL	10
TECH	INICAL SPECIFICATION	11
OVEF	NIEW	12 -14
INST	ALLATION	15-17
OPEF	RATION	18-20
CLEA	NING & CARE	21-22
ROUT		23-24
PART	'S LIST	25
TROL	JBLESHOOTING	26-29
ЕМС	STATEMENT	30-34
GUAF	RANTEES & WARRANTIES	35
DECL	ARATION OF CONFORMITY	36-37

WARNINGS & CAUTIONS



READ THIS INSTRUCTION MANUAL AND OBSERVE SAFETY INSTRUCTIONS.



WARNING

- Disconnect the battery when not in use. Replace the battery if less than 12 Volts. Work should only be undertaken by an authorised person.
- This system must be properly installed and operated as directed by this user manual.
- The system should be checked regularly to ensure correct operation. Loss of function will remove all pressure relieving properties that this system provides.
- This system is intended for use as part of a pressure ulcer prevention program; do not rely solely on this device to achieve the result. The medical professional is responsible for applying best medical judgment when using this system.
- Select the correct setting for the occupant's weight and therapy required (see pages 18-19). Care should be taken not to accidently change pressures once set as the effectiveness of the therapy may be reduced.
- In order for alternating air pressure range to be effective, avoid placing objects on the surface that may obstruct the movement of air between the cells. For the same reason, discourage people from sitting on the edge or on the end of the mattress whilst it is in use.
- All hoses must be free of kinks, twists and must be properly connected and positioned so as not to cause any obstruction.
- Do not position the system in a way that prevents access to the disconnection device (mains power plug).
- Ensure the mains lead or pump cannot become trapped or crushed, e.g. by raising or lowering of bed or bed rails or any other moving object.
- Check the mains lead is damage free and positioned so as not to cause an obstruction, or injury, e.g. Strangulation or trip hazard.
- Ensure that the electricity supply is of the type stated on the pump unit.
- Protect your system from open flames. Ensure that the system is not used in the presence of flammable anaesthetics.
- Do not place device on or near a heat source or cover pump with bedding.
- Harvest Healthcare advise against smoking whilst the system is in use, to prevent the accidental ignition of associated items which may be flammable, such as bed linen.

WARNINGS & CAUTIONS

- Do no expose the pump to liquids.
- Do not use with hot water bottles or electric blankets.
- Wireless equipment such as mobile phones should be kept at least 10ft / 3m away from the system.
- Do not allow sharp objects to puncture the mattress material.
- The mattress and pump should be cleaned between patient uses.
- Do not use bleach, chlorine releasing agents in concentrations over 1000 ppm, solvents or alcohol-based cleansers, e.g. Phenicol, Hibiscrub, Clearsol, Stericol, Hycoline as these will destroy the mattress material. Full cleaning instructions can be found on **page 21-22**.
- Suitable for continuous use.
- Do not modify the mattress or pump unit in any way.
- Do not connect to any other medical device or equipment.
- Not for use in an oxygen enriched environment.
- Not for use in an outdoor environment.
- Store the system in a clean and dry environment, out of direct sunlight.



Electrical equipment can be hazardous. Only authorised technical personnel should remove the rear pump case for maintenance. Removal of the case by unqualified personnel will invalidate the warranty.



Before cleaning the unit ensure that the electrical supply to the pump has been disconnected by removing the plug from the power supply.



Do not use this system for lifting the patient. This will damage the system and could put the patient at risk.



This product is fire rated. The mattress cover material on the mattress is tested to BS7175:1989 Crib 5. The internal cells are tested to BS EN 597-1:1995. Use of this product should be subject to a risk assessment in which all hazards are considered.

BEFORE USING THIS SYSTEM FOR THE FIRST TIME:

- Read through this instruction manual conscientiously from start to finish.
- Please note that the various safety instructions must be observed.

Harvest Healthcare products bear the CE mark and meet all safety and functionality requirements.

These safety requirements can only be met if the user is satisfied with the proper condition of the product (including accessories) before use.

GENERAL INFORMATION

The **Royal System** is an alternating pressure relieving mattress system used in the prevention and treatment of pressure ulcers, and is recommended for use by a patient who is at risk from developing pressure sores. The mattress is fitted with a vapour permeable two way stretch cover.

By using the established principles of alternating therapy, the **Royal System** offers the patient comfortable and relaxing support that can both prevent tissue breakdown and enhance healing.

The **Royal Pump** unit is lightweight and compact; its features include an audible and visual low-pressure warning, intergrated battery backup and a manual pressure / comfort control function.

The **Royal Mattress** is made up of 17 individual cells. 14 alternating air cells and 3 static cells at the head end. The alternating cells are split into 2 sections - odd cells e.g. 1,3,5 etc and even cells e.g. 2,4,6 etc. These two sections will alternate through a 10 minute cycle in which time both sets of air cells will have inflated and deflated sequentially. All air cells are individually replaceable should any damage occur.

The quick release 3 pipe connector complete with transport cap enables easy patient transport arrangements. For rapid deflation of the system simply twist open the CPR.

The **Royal Extra Mattress** is based on the Royal Replacement Mattress System with an additional foam surround, suitable for wider beds The mattress has 14 alternating air cells and 2 static cells at the head end. Available in 1200mm and 1370mm width, the mattresses have an increased maximum weight limit of 50 stones / 317 kg. The surround is visco topped for comfort and can be used on a profiling bed.

1 DEFINITION OF THE GROUPS MENTIONED

OPERATOR

An operator is any natural or legal person who uses the equipment or on whose instruction it is used (e.g. nursing homes, specialised retailers, health insurance companies, medical suppliers).

USER / CARE PERSONNEL

Users are persons who as a result of their vocational training, experience or briefing are authorised to operate the equipment.

Furthermore, the user/ care personnel can recognise and avoid potential dangers and assess the clinical condition of the service user.

PATIENT / OCCUPANT / SERVICE USER

The person in need of care, handicapped or infirm.

QUALIFIED PERSONNEL

Qualified personnel are employees of the operator who as a result of their vocational training or briefing are entitled to deliver, assemble, disassemble and transport the product.

2 NON-COMPLIANT USE

All uses deviating from the intended purpose, which may also be hazardous as a result. This includes for example:

- Incorrect installation.
- Operation by persons who have not been instructed in its use.
- Using the system with non-approved parts/accessories.
- Using the system if any of the components are damaged or faulty.

3 SAFETY INSTRUCTIONS

3.1 GENERAL SAFETY INSTRUCTIONS



During the briefing, specific attention must be drawn to any potential dangers which can occur despite correct operation. Before putting the product into service for the first time, the Instruction Manual must be read conscientiously and in detail by the user / care personnel.

Programming of the system should be carried out by competent trained personnel. Use only spares and accessories supplied or approved by Harvest Healthcare. Only suitably trained personnel are allowed to operate the system.



The mains cable must be free and not be allowed to be caught up in the bed's moving mechanisms. The mains cable may be torn out of its strain relief and damaged or it may be pulled out of its socket and electric leads exposed as a result.

If the mains cable or the mains plug are damaged, the relevant part must be replaced. This work should be carried out by the manufacturer or authorised service agents.

When connecting the mains plug do not use multiple sockets since liquids may penetrate into these (fire hazard and electric shock).

3.2 SAFETY INFORMATION FOR THE OPERATOR

With the help of this Instruction Manual, instruct each user in the safe operation of this system before it is put into service for the first time.

Advise the user of any hazards which may occur if not handled correctly.

Only persons who have been properly instructed may operate this system. This also applies for persons who only operate the system on a temporary basis.

3.3 SAFETY INFORMATION FOR THE USER

Ensure that the operator instructs you in the safe operation of this system.

In addition, pay particular attention to the Warnings and Cautions (**page 4-5**) and the general safety information as described in **3.1**.

If there is a suspected fault or damage, unplug the mains plug from the socket and follow the power down procedure (**page 19**). Clearly mark "Out of Order" and take out of service immediately, and inform the person in charge without delay.

3.4 SYMBOLS USED

This symbol indicates general hazards. There is a danger to life and health. Danger acid, including battery acid. Danger risk of explosion, Conformity mark in accordance with the Medical Device Directive (93/42 EEC). The symbol for Protection Class II device, double insulated. The symbol for type B device according to EN 60601-1. Handle with care This way up Keep dry Recycling symbol. Refers to packaging that can be recycled (cardboard) Fragile, handle with care This product must be disposed of in a designated refuse bin for waste electronic devices (WEEE) in the European Union. Do not dispose of as normal domestic waste. Symbol for direct current. This product may only be used indoors. IP: Ingress protection (protection from vertically dripping water) Wear protective goggles. No smoking. No naked flames.

Children should not be allowed to play with or attempt to open this battery.

Read instructions / consult manufacturers guide

3.5 CLEANING & DISINFECTION

Do not immerse electrical components in water but wipe with a damp cloth only. The electrical components must not be cleaned with a high-pressure cleaner or water jet. Disinfection by wiping only is allowed.

Full cleaning and disinfection instructions can be found on **pages 21-22**.

3.6 SERVICING & MAINTENANCE

Disconnect the battery when not in use. Replace the battery if less than 12 Volts. Work should only be undertaken by an authorised person.

Servicing must only be carried out by qualified personnel.

A technical check and/or safety inspection must be conducted at least once a year and after a lengthy break in use.

Any defects, damage or signs of wear must be rectified without delay. Only original spare parts from Harvest Healthcare may be used, otherwise all guarantees or warranties will be excluded.

See pages 23-24 for Routine Maintenance and Servicing.

3.7 SERVICE LIFE & DISPOSAL



The system must not be disposed of as normal domestic waste after its service life, but must be disposed of in a designated refuse bin for waste electronic devices (WEEE) in the European Union. Do not dispose of as normal domestic waste.

Our Full Terms & Conditions including product warranties are available by request or can be found on our website **www.harvesthealthcare.co.uk.**



PARTS AND DATA MAY UNDERGO FURTHER DEVELOPMENT AND THEREFORE DEVIATE FROM THE DETAILS GIVEN.

TECHNICAL SPECIFICATION

ROYAL SPECIFICATION

Product Code Pressure Sore Risk Level Minimum Patient Weight Maximum Patient Weight Inflated Mattress Dimensions Mattress Weight Fire Retardancy (Cover) HR400 Very High Risk 5 Stone / 31.7 kg 39 Stone / 247 kg 1940 x 920 x 205 mm 9 kg BS7175 Crib 5

ROYAL EXTRA SPECIFICATION

Product Code Pressure Sore Risk Level Minimum Patient Weight Maximum Patient Weight Inflated Mattress Dimensions Fire Retardancy (Cover) **1200** - HR430C / **1370** - HR430D Very High Risk 5 Stone / 31.7 kg 50 Stone / 317 kg 1940 x 1200 / 1370 x 240 mm BS7175 Crib 5

PUMP TECHNICAL SPECIFICATION

Pump Model No. NP-102 Pump Dimensions 240 x 370 x 120 mm Pump Weight 5.2 ka **Operating Cycle** 10 minutes 0.194A / 44.5VA Power Rating Voltage 230V / 50Hz Noise Level NC30 Fuse T 1AH 250V Medical Classification Type B Applied Part IP Rating IPX1 EN 60601-1. EN 60601-1-2 Safetv Standards **Electric Shock Protection** Class II ME equipment externally powered/ internally powered.

OVERVIEW

ROYAL OVERVIEW

CPR Valve

The mattress can be rapidly deflated through the use of the CPR valve, allowing emergency personnel to begin resuscitation.

Royal Cells

The Royal has 17 independent cells which can be individually removed and replaced to allow for cost effective repair and in-depth cleaning.

Mattress Cover

The cover comprises of a two-way vapour permeable and fully welded stretch PU top with a durable base fabric.

Royal Pump

The Royal pump has a moulded ABS case with non-slip feet on the base and integrated bed hooks.

Feed Tubes The Feed tubes are flexible, durable and have excellent anti-kink properties.

SHOULD HAVE

OVERVIEW

ROYAL EXTRA OVERVIEW

CPR Valve

The mattress can be rapidly deflated through the use of the CPR valve, allowing emergency personnel to begin resuscitation.

Royal Cells

The Royal Extra has 16 independent cells which can be individually removed and replaced to allow for cost effective repair and in-depth cleaning.

Mattress Cover

The cover comprises of a two-way vapour permeable and fully welded stretch PU top with a durable base fabric.

Foam Surround

The Royal Extra is surrounded by a foam wall. This makes the Royal Extra suitable for wider beds and a higher maximum weight.

Feed Tubes

The Feed tubes are flexible, durable and have excellent anti-kink properties.

Royal Pump

The Royal pump has a moulded ABS case with non-slip feet on the base and integrated bed hooks.

OVALEXTRO



PUMP OVERVIEW





INSTALLATION

INSTALLING THE ROYAL SYSTEM

- **1** Remove the mattress from its packaging and lay the parts out on the floor. You should have the following items:
 - Carry bag (holdall)
 - Mattress with feed tubes attached
 - Pump
 - Instruction booklet




If you intend to keep this mattress system in storage at some point please retain the packaging. This will lengthen the life of the mattress.



Prior to installing the mattress, check that there are no protruding/ sharp objects which may puncture the cover or air cells.


- 2 Carefully unroll the mattress on to the mattress platform. Ensure that the pipes at both the head and foot ends are kink free and straight to prevent restriction of air flow. Ensure that the air tubing to the pump is at the foot of the bed.
- 3 The straps on the under side of the mattress are used to secure the mattress to the bed. Secure the straps using the plastic clips around the bed frame and tighten.
 Note: Secure the straps to the moving parts of the platform only.

Care must be taken when fixing the mattress to the bed frame.



INSTALLATION

6 The pump has integral bed hooks for hanging the pump on the foot end of the bed.



7 Check that the CPR valve is set to the **CLOSED** setting as indicated below.





The mattress will not inflate if the CPR valve is open.

8 Connect the feed pipes to the pump using the quick release coupling and ensure the connection has securely clicked into place.



9 When the mattress is ready to be inflated, insert the mains plug into the wall socket and turn on the power.

INSTALLATION

INSTALLING THE ROYAL EXTRA SYSTEM (ADDITIONAL INFORMATION)

Please refer to the previous installation instructions (**pages 15-16**). However, please be aware of the following additional guidelines when installing the Royal Extra systems:

- Your Royal Extra system will be delivered to you in a box. This prevents the foam surround from damage during trasportation. In addition, any attempt to roll or fold the mattress may result in damage to the foam surround
- The Royal Extra mattress does not include securing straps. Please ensure that the mattress is correctly sized for the bed and that the mattress is positioned onto the mattress base using the mattress guides.
- After switching the system on, It is advised that the system is left to fully inflate before
 putting into use.



If you intend to keep this mattress system in storage at some point please retain the packaging. This will lengthen the life of the mattress.



Ensure that the Royal Extra Mattress is fully inflated before being put into use.

OPERATION

INFLATING THE MATTRESS

1 Switch on: The pump will automatically begin to inflate and start a 30 minute mattress inflation mode until the correct pressure is reached.

The LCD display will read **Please Wait** until the correct pressure is acheived.

2 When the mattress has fully inflated the LCD display will read **Mattress Ready**. The Service User can now be layed on the mattress; this will automatically begin an Auto Setting function, in which the pump will adjust its pressure settings to suit the user's requirements.

Alternitavely, you can manually set the comfort setting (please see comfort control gude on the next page).



Set the pump to run at the correct pressure to suit the weight of the service user. Refer to the Comfort Control Guide (see next page) or on the front of the pump.

IN THE EVENT OF A POWER FAILURE

- 1 Should the power supply to the pump be cut, the system will automatically switch to battery back-up.
- 2 If the power is restored, the pump will continue on the selected setting.
- **3** The battery will opperate for approx 6-8 hours, afterwhich the system will deactivate and will therefore need to be connected again to a live mains power supply.

DEFLATING THE MATTRESS / CPR VALVE



If rapid deflation is required, simply twist the CPR to the $\underline{\text{OPEN}}$ position.

To deflate the mattress, simply open the CPR valve and disconnect the pump from the mattress using the quick release connector.



OPERATION

SWITCHING OFF THE SYSTEM

- 1 Turn off and unplug from the mains supply (the battery backup will operate).
- **2** Press the function button and hold continuously until the unit has turned off. (approx. 30 seconds).

TRANSPORT MODE

A transport cap is provided to hold the air in the mattress if the pump is disconnected temporarily for any reason.

To prevent deflation please ensure that the transport cap (attached to the air tubes) is fitted to the 3 pipe connector on the mattress (see picture)






COMFORT CONTROL GUIDE

	WEIGHT KG	WEIGHT LB	WEIGHT STONES LB
LOW	32 - 75 kg	70.5 - 167 lb	5 - 11 stone 11 lb
MEDIUM	76 - 100 kg	166 - 220 lb	12 - 15 stone 10 lb
HIGH	100 - 254 kg	220 - 560 lb	15 stone 10 lb - 40 stone

OPERATION

OPERATION OF THE SYSTEM



RED LOW PRESSURE WARNING LIGHT

The low pressure light will glow red when a fault is detected. Check that the CPR is closed and ensure that all pipes are connected to the mattress.

If the warning light does not go out, contact your supplier for assistance.

ADJUSTING THE WEIGHT SETTING MANUALLY

To change the weight setting manually, press and hold the function button for 12 seconds, then release. The next weight setting will then engage. Repeat this operation until you have reached the required weight setting (i.e. light, medium or heavy).

ADJUSTING THE COMFORT SETTING

Adjust the comfort setting by pressing the function button to cycle the pump through the different settings i.e. Medium (default setting turning the system on) Static, Low, Medium and then High.

STATIC FUNCTION

This is used during patient care and inflates all 30 cells for 20 minutes; after this the pump will return to its previous setting.

CLEANING & CARE

WARNING

Ensure that the mains power supply to the pump is disconnected before cleaning

Eye protection, gloves and protective clothing should be worn when carrying out cleaning and disinfection procedures

When disinfecting the system, Harvest Healthcare recommends the following guidelines which have been developed to comply with recognised infection control procedures. These procedures are also to be used to prevent cross infection when transferring the system between patients.

MATTRESS

During general use the mattress and internal tubes can be cleaned by wiping with a mild detergent solution.

Where necessary the mattress cover can be removed for laundering or sterilisation. Where there is staining or body fluids on the mattress, cells or tubing, wash thoroughly with soap and water, then wipe with a sodium hypochlorite solution diluted to 1000ppm before laundering.

Mattress covers may be laundered as follows:

- 1 Pre-wash Cold 10 minutes
- 2 Main Wash 80°C 10 minutes
- **3** Followed by cold rinses and extraction.



Do not use abrasive cleaners, phenol disinfectants, solvents or alcohol-based cleansers, e.g. Dettol, Phenicol, Hibiscrub, Clearsol, Stericol, Hycoline, as these will damage the cover materials.

Do not iron.

Ensure that the mattress/cushion is thoroughly dried before remaking the bed or placing in storage.



HYPERCARBONATE AND PHENOL-BASED SOLUTIONS SHOULD NOT BE USED AS THESE WILL DAMAGE THE MATTRESS COVER

CLEANING & CARE

PUMPS

For general cleaning wipe with a soft cloth dampened with a mild detergent and water solution. This may be followed by either wiping with a sodium hypochlorite solution to a dilution of 1000ppm (parts per million) or by using alcohol wipes.



DO NOT USE HYPERCARBONATE, PHENOL-BASED CLEANING SOLUTIONS, ABRASIVE COMPOUNDS OR CLEANING PADS.

NOTES

- Following the use of a detergent and or disinfectant solution rinse the mattress cover with clean water using a clean cloth and allow to dry.
- Frequent or prolonged exposure to high concentrations of disinfectant solutions will reduce the useful life of the mattress cover.
- Where high concentration disinfectants e.g. > 10,000ppm chlorine releasing agent (e.g. Haztab or bleach) or combined cleaning/chlorine releasing agent (e.g. Chlorcleam, Actichlor) and detergent solutions are used to remove blood or other body fluids, mattresses should be thoroughly rinsed with clean water to remove any residues. This will help prevent any long term compatibility issues associated with disinfectant residues.
- Alternatively, disinfection may be achieved by laundering at temperatures not exceeding 80°C for 10 minutes which may include a chlorine rinse.

TRANSPORT & STORAGE

Storage conditions as follows:

-15 °C without relative humidity control; and +40 °C at a relative humidity up to 93%, non-condensing. An atmospheric pressure range of 700 hPa to 1 060 hPa. Suitable for all standard modes of transport when in the correct packaging.

Operation Conditions:

A temperature range of +5 °C to +35 °C; A relative humidity range of 15% to 93%, noncondensing; and Operational Atmospheric Pressure: 700 hPa to 1060 hPa. Suitable for pollution degree 2. Operational altitude \leq 2 000 m.

Transportation of the mattress system:

The mattress should be loosely rolled lengthwise with the cover innermost, taking care not to strain the feed pipes. It can then be stored / transported in the carry bag with the pump, mains cable and this booklet. Do not stack bagged mattresses more than two high.

ROUTINE MAINTENANCE

These checks should be carried out at each decontamination process, i.e. between patients or patient occupancy and weekly for longer term patients.

MATTRESS

The mattress cover, which is made from waterproof and vapour permeable material, should be kept clean. Take care to avoid puncturing cover with sharp objects whilst performing the maintenance checks:

- 1 Remove cover and inspect for damage, tears or staining, which could lead to contamination of the internal parts.
- 2 Check that the zips are sound and in good working order.
- 3 Check that all connectors are fitted properly to prevent leaking of air.
- 4 Check that all cells are attached to the base sheet by the pop fittings provided.
- 5 Check the stitching on the straps and the seams to ensure no tearing or fraying has occurred.

PUMP

- 1 Check the pump casing for cracks or other damage that could be dangerous.
- 2 Check the power cord (ensure there are no bare wires).

If any faults are detected report to your distributor for replacements to facilitate repairs.

COMPONENTS

- Check air cells and mattress interior for signs of damage or contamination, e.g. staining or fluid ingress at each decontamination process, i.e. between patients or patient occupancy (or weekly for longer term patients).
- The individual cells can be wiped clean with a mild antiseptic solution.
- All cells are replaceable and can be sourced from Harvest Healthcare.

POWER UNIT

Disconnect the power unit from the electricity supply before carrying out maintenance, repairs or cleaning.

Check all electrical connections and power lead for signs of wear and damage. The power unit can be wiped down with detergent, disinfectant solution or wipe*. Do not use solvents. Unsuitable for sterilisation.

ROUTINE MAINTENANCE

- * In line with the MHRA Medical Device Alert (MDA/2013/019), Harvest Healthcare advises customers to use pH neutral, high-level disinfectant cleaning products to sanitise reusable medical devices to prevent damage to materials and the degradation of plastic surfaces after prolonged use.
- The use of inappropriate cleaning and detergent materials on medical equipment could damage surfaces and may compromise the ability to decontaminate medical devices adequately or may interfere with device function

At end of use dispose of the pump / mattress in accordance with the local regulations including WEEE requirements, which apply to the pump and SMPS only.

SERVICING YOUR SYSTEM

The Royal system should be serviced every **12 months** by Harvest Healthcare approved personnel using genuine Harvest Healthcare spare parts.



Failure to follow the Royal service schedule may invalidate future warranty claims (Gurantees & Warranties can be found on page 35).

PARTS LIST

APPLIED PARTS

HR402	Royal Mattress (only)	
NP-102	Royal Pump (only)	
HR423	Royal Top Cover and Base	
HR433C	Royal Extra 1200 Top Cover and Base	
HR433D	Royal Extra 1370 Top Cover and Base	
HR432C	Royal Extra 1200 Foam Surround	
HR432D	Royal Extra 1370 Foam Surround	
HS384	Royal Inner (Complete)	
HP574/3	CPR Valve	
ННЗЗ	Seat Cushion (Optional)	

REPLACEMENT PARTS

HS382	Royal Cell (Static)	
HS383	Royal Cell (Alternating)	
HS381	Inner Base Sheet	
HP468_1	3 Pipe Connector (only)	
HP611	10mm T Connector	
HP589	10mm Elbow	
HLDW300	Pump Hangers	
HP468	Transport Cap with Tag	
HR405B	Royal Bag/ Holdall	

TROUBLE SHOOTING

FAULT	СНЕСК ТНАТ	STAGE 2 CHECK	IF PROBLEM PERSISTS
Pump shows no indication that it is powered up	² ump shows no indication 1. Mains plug is plugged in and hat it is powered up power switched on.	 Connect the pump to the nearest (working) mains outlet. 	Contact Harvest Healthcare technical support.
	 The power switch on the pump is switched on. 	Replace the Fuses with the correct T1A fuses.	Before calling:
	 The fuse in the mains plug is not blown. 	 Try a different device in the mains outlet 	Please ensure you have the serial number and model of equipment.
	 The wall socket that the pump is connected to is working correctly 		Please record details of the results of the recommended tests. (Notes pages are provided at the back of this user manual).

FAULT	СНЕСК ТНАТ	STAGE 2 CHECK	IF PROBLEM PERSISTS
The Pump appears to be running but the mattress is not inflating correctly and or the low pressure light is illuminated.	PLEASE NOTE Inflation can take up to 30 minutes. If the Low Pressure Alarm sounds, press the FUNCTION button.	 Disconnect and then re- connect the hoses to the outlet on the side of the pump. Open then reclose the CPR valve, make sure the valve is not trapped in the bed 	Contact Harvest Healthcare technical support. Before calling: Please ensure you have the serial number and model of equipment.
	 The hoses are routed correctly (not kinked) and connected to the pump correctly. The CPR valve is not trapped and is in the closed position. 		Please record details of the results of the recommended tests. (Notes pages are provided at the back of this user manual).
	 There are no leaks in the mattress. The tubes in the mattress are not disconnected or kinked. 	 Replace any damaged or leaking mattress parts with the correct genuine Harvest Healthcare spare parts. Straighten out any kinked pipes and reconnect any disconnected joints. 	

TROUBLE SHOOTING

FAULT	СНЕСК ТНАТ	STAGE 2 CHECK	IF PROBLEM PERSISTS
Some of the cells appear to	This is normal for alternating pressure therapy.	sure therapy.	
טם ממומנים.	The mattress is made up of indivic cells e.g. 1,3,5 etc and even cells which time both sets of atternating	The mattress is made up of individual air cells. The alternating section is split up into 2 cells e.g. 1,3,5 etc and even cells e.g. 2,4,6 etc. These two sections will alternate thro which time both sets of alternating air cells will have inflated and deflated sequentially.	The mattress is made up of individual air cells. The alternating section is split up into 2 sections consisting of odd cells e.g. 1,3,5 etc and even cells e.g. 2,4,6 etc. These two sections will alternate through a 10 minute cycle in which time both sets of alternating air cells will have inflated and deflated sequentially.
The system does not appear to be alternating.	 Carefully mark one of the inflated cells with a pen. 	 Monitor the cell for 7 minutes to see if it deflates. 	Contact Harvest Healthcare technical support.
	2. Ensure that there are no kinks in the ninework down	2. Straighten out any kinked	Before calling:
		3. Press the function button until	Please ensure you have the serial number and model of equipment.
	 Ensure that the static function on the pump has not been activated 	the static light goes out and set the pressure control to the annicoriate setting for the	Please record details of the results of the recommended tests. (Notes pages are
		person on the mattress.	provided at the back of this user manual).

FAULT	снеск тнат	STAGE 2 CHECK	IF PROBLEM PERSISTS
The pump is vibrating or making excessive noise.	The pump is fitted to the bed correctly	Reposition the pump unit.	Contact Harvest Healthcare technical support.
The mattress is uncomfortable.	Check the comfort setting on the pump.	Set the pump to the correct setting using the guide on the front of the pump case.	Before calling: Please ensure you have the serial number and model of equipment. Please record details of the results of the recommended tests. (Notes pages are provided at the back of this user manual).

A comprehensive service manual for this system is available upon request.





USE OF PORTABLE TELEPHONES OR OTHER RADIO FREQUENCY (RF) EMITTING EQUIPMENT NEAR THE SYSTEM MAY CAUSE UNEXPECTED OR ADVERSE OPERATION

The Royal has been tested to EN60601-1-2:2007, regarding its ability to operate in an environment containing other electrical/electronic equipment (including other medical devices).

The purpose of this testing is to ensure that the Royal is not likely to adversely affect the normal operation of other such equipment and that other such equipment is not likely to adversely affect the normal operation of the Royal

Despite the testing of the Royal that has been undertaken, normal operation of the Royal can be affected by other electrical/electronic equipment and portable and mobile RF communications equipment.

As the Royal is medical equipment, special precautions are needed regarding EMC (electromagnetic compatibility).

It is important that the Royal is configured and installed/put into service, in accordance with the instructions/guidance provided herein and is used only in the configuration as supplied.

Changes or modifications to the Royal may result in increased emissions or decreased immunity of the harvest Royal in relation to EMC performance.

The Royal should not be used adjacent to or stacked with other equipment. If adjacent or stacked use with other equipment is necessary, the Royal and the other equipment should be observed / monitored, to verify normal operation in the configuration in which it will be used.

For the purposes of EN60601-1-2 the Royal has an essential performance (the Royal must maintain required pressure within a tolerance of +/-10%).

- 1 The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of the equipment or system as replacements for internal components, may result in increased emissions or decreased immunity of the equipment or system.
- 2 The use of the accessory, transducer or cable with equipment and systems other than the Royal may result in increased emissions or decreased immunity of the equipment or system.

EMC STATEMENT

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

The Royal is intended for use in the electromagnetic environment specified below. The customer or the user of the harvest Royal should assure that it is used in such an environment.

EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE	
RF emissions CISPR 11	Group 1	Uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR 11	Class B	Suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that	
Harmonic emissions IEC61000-3-2	Class A	low-voltage power supply network that supplies buildings used for domestic purposes.	
Voltage fluctuations / flicker emissions IEC61000-3-3	Complies		

EMC STATEMENT

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

The Royal is intended for use in the electromagnetic environment specified below. The customer or the user of the harvest Royal should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Electrostatic discharge (ESD) IEC61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	\pm 1 kV line(s) to line(s) \pm 2 kV line(s) to earth	\pm 1 kV line(s) to line(s) \pm 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11			Mains power quality should be that of a typical commercial or hospital environment. If the user of the harvest Royal requires continued operation during power mains interruptions, it is recommended that the harvest Royal be powered from an uninterruptable power supply or a battery.
Power frequency (50/60Hz) Magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.

NOTE: U_{τ} is the a.c. mains voltage prior to application of the test level.

EMC STATEMENT

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

The Royal is intended for use in the electromagnetic environment specified below. The customer or the user of the harvest Royal should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE	
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Royal including any cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.	
Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	Recommended separation distance (d)	
			d = 1.2√P	
			d = 1.2√P 80 MHz to 800 MHz	
			d = 2.3√P 800 MHz to 2.5 GHz	
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).	
			Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range.b	
			Interference may occur in the vicinity of equipment marked with this symbol:	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the harvest Royal is used exceeds the applicable RF compliance level above, the harvest Royal should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the harvest Royal

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.



GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

The Royal is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Royal can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the harvest Royal as recommended below, according to the maximum output power of the communications equipment.

RATED MAXIMUM OUTPUT POWER OF TRANSMITTER	SEPARATION D	ISTANCE ACCORDIN OF TRANSMITTER (
(W)	150 KHZ TO 80 MHZ D = 1.2√P	80 MHZ TO 800 MHZ D = 1.2√P	800 MHZ TO 2.5 GHZ D = 2.3√P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures objects and people.

GUARANTEES & WARRANTIES

MATTRESS (COVER AND INTERIOR COMPONENTS)

All Harvest Healthcare Ltd Mattress are covered by warranty for a period of 12 months from date of purchase. Damage through incorrect use and penetration by sharp instruments will invalidate this warranty.

PUMP

The Royal Pump is covered by warranty for a period of 3 years from the date of purchase. This excludes all serviceable parts such as the bellows and filters which are recommended to be changed every 12 months in line with the service schedule.

A service manual for this system is available upon request.

GUARANTEE

Harvest Healthcare Ltd guarantees to repair or replace all goods supplied to its customers which are found to be defective whilst still in their applicable warranty period. All warranties are subject to the following conditions:

- **a** Warranty/ guarantee is subject to all guidelines being adhered to.
- **b** That the equipment has been used for the purpose for which it was intended.
- **c** That the usage has been on a fair wear and tear basis. This does not include user damage.
- **d** That Harvest Healthcare Ltd's cleaning/ disinfecting guidelines have been followed.
- e Harvest Healthcare Ltd's maintenance guidelines have been followed (Please refer to the product manual).
- **f** That ALL maintenance has been carried out by a suitably qualified and competent person.
- **g** That all parts used are OEM (Original Equipment Manufacturer) parts and were supplied by Harvest Healthcare Ltd either directly or through a distributor.
- **h** All warranties begin from the time the product leaves the premises of Harvest Healthcare Ltd.
- i All repairs and replacements will be at the sole discretion of Harvest Healthcare Ltd.

Our standard terms and conditions of sale can be found on our website or by request to Harvest Healthcare Ltd

DECLARARTION OF CONFORMITY

Declaration of Conformity Annex VII EU Directive 93/42/EEC

We, as company:

Harvest Healthcare Limited Sheaf House Bradmarsh Way Bradmarsh Business Park Rotherham, S60 1BW

Confirm on our own behalf that the medical product: ROYAL ACTIVE MATTRESS SYSTEM

Complies with all applicable requirements in Appendix I of the EU directive 93/42/EEC.

The following compliance evaluation process was applied: Annex VI

In the event of modification of this product without consultation with the manufacturer, this declaration of conformity will lose its validity.

Rotherham, 01.12.2015

Dree Mluna

Director

DECLARARTION OF CONFORMITY

Declaration of Conformity Annex VII EU Directive 93/42/EEC

We, as company: Harvest Healthcare Limited Sheaf House Bradmarsh Way Bradmarsh Business Park Rotherham, S60 1BW

Confirm on our own behalf that the medical product:

ROYAL EXTRA ACTIVE MATTRESS SYSTEM

Complies with all applicable requirements in Appendix I of the EU directive 93/42/EEC.

The following compliance evaluation process was applied: Annex VII

In the event of modification of this product without consultation with the manufacturer, this declaration of conformity will lose its validity.

Rotherham, 01.12.2015

me Mh

.....

Director

NOTES

NOTES



harvesthealthcare[®]

Sheaf House, Bradmarsh Way, Bradmarsh Business Park, Rotherham, S60 1BW

 T +44 (0)1709 377 172
 F +44 (0)1709 377 173
 E sales@harvesthealthcare.co.uk www.harvesthealthcare.co.uk

Serial No:

DOCUMENT REFERENCE: HR400/MANUAL02 - May 2017

The Royal pump supplied by Harvest Healthcare Ltd. complies to the standard EN 60601-1.

© Copyright Harvest Healthcare Harvest Healthcare Limited. Company No: 07210261 Our Full Terms & Conditions are available by request or can be found on our Website Harvest Healthcare reserves the right to alter or amend this document without notice