

WOBURN ULTIMATE PROFILING BED



Product Code:

HLB798/12

Description:

The Woburn Ultimate Profiling Bed is our higher capacity bed within the Woburn range. Featuring triple aluminium side rails and under-bed strengthening, the Ultimate excels in patient support and ease of infection control. Available in 1200mm wide, with a 60 Stone weight capacity, this bed is able to support most bariatric patients. Comfort adjustments are made with ease due to the powerful lifting components and the high quality electronic control system.

Safe Working Load:

380kg (60 Stone)

Max. Patient Weight:

300kg

Dimensions:

2195 x 1340mm (1200mm mattress base)

Max. Mattress Depth:

260mm

Adjustable Height:

275 - 690mm

Angle of Back Rest:

Up to 70 degrees

Angle of Knee Break:

Up to 30 degrees

Bed Weight:

157kg + 16kg with rails

Side Rails

Description:

Aluminium side rails with plastic end caps

Dimensions:

1973 x 70 x 28mm

Key Information:

4 Section Profiling
Electronically Powered Back Rest & Thigh Rest
Manual Leg Rest Adjustment
Wraparound Wooden Head and Foot Boards
Aluminium Side Rails
Auto-Regression as standard
Trendelenburg Function as standard
200mm Extension - Optional
Lockable Handset
Low Voltage – 29v; 9v Battery Backup
Comprehensive Manual
CE marked

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

Important Safety Notice

<Ufj Ygh< YUth.VWYF @X gi dd]Yg Yei]da YbhUbX Udd]UbWg hc Wj Yf U fub[Y cZa YX]WU' VebX]jicbg" GdYVW]Wk.c]M' cZYei]da Ybhig h.Y gc.Y fYgdcbg]V]]mcZh.Y Wb]WU' dfUW]jicbYf cf WUWYF
cj YfgYY]b[]bX]]]X] U' dU]YbhWUWY" Rig_ UggYgga Ybhg gdYVW]Wk' i gYfa i ghVY WUff]X ci hd]cfhc h.Y i gY'cZU" dfcXi Wtg"
PrcXi Wtg UobX XUHU Veb]bi U'mi bXYf[c XYj Ycda YbhUbX a Umh.YfYz:fy XYj]U]Y Zca h.Y XYHU]g [] Yb"

Version 2: Jan 2018