

# HARVEST 1 PUMP



**Product Code:** HPU1

**Description:**

- Harvest 1 is our manually operated pump, ideal for use in community and care home settings
- Simple interface with a dial
- Easy to connect and disconnect and is lightweight, creating ease of use for the carer
- Has both visual and audible alarms for safety
- Fits with the Kensington overlay, Sandringham Replacement mattress Hampton Replacement mattress and Prime Comfort Active Mattress

**Dimensions:** 90 x 150 x 290mm

**Weight:** 1.7kg

**Electrical Safety Standards:** EN 60601-1: 2006 and A11:2011

**Voltage:** 230V 50Hz

**Pressure Range:** 30-80 mmHg

**Variable Cycle Time:** 10 minutes

**Fuse:** T1A 250v

**Warranty**

In terms of Harvest Healthcare Guarantee and Warranties the Harvest 1 Pump is guaranteed for one year from date of purchase.

**Supplier Contact Details:**

t: +44(0)1709 377172 e: [enquiries@harvesthealthcare.co.uk](mailto:enquiries@harvesthealthcare.co.uk) w: [www.harvesthealthcare.co.uk](http://www.harvesthealthcare.co.uk)  
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

**Important Safety Notice**

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

**Version 2: March 2018**