

HARVEST 3 PUMP



Product Code:

HPU3

Description:

- Harvest 3 is our high-end digital pump, ideal for use in all settings including acute
- Caters for bariatric products and more complicated clinical needs
- Easy to use
- L.E.D display, lock out and static mode
- Maxi-firm mode and four cycle times
- Visual and audible alarms for low pressure and power failure
- Very high capacity, fitting all of our active mattresses including our bariatric range

Dimensions:

340 x 125 x 220mm

Weight:

3.5kg

Electrical Safety Standards:

EN 60601-1: 2006 and A11:2011

Voltage:

230V 50Hz

Pressure Range:

20-65 mmHg

Variable Cycle Time:

10 / 15 / 20 / 25 minutes

Fuse:

T1A 250v

Warranty

In terms of Harvest Healthcare Guarantee and Warranties the Harvest 3 Pump is guaranteed for three years from date of purchase.

Exclusions:

1. Components of pumps which are recommended for annual service, are guaranteed for one year.

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 2: March 2018