




harvesthealthcare®
equipped to care

Static Cushions



www.harvesthealthcare.co.uk

POWER TO
DELIVER

HARVEST HEALTHCARE
YOUR FIRST CHOICE...



TOTAL SUPPORT TEAM

We offer a full package to support you from installation to service and repair.



HARVEST FAST TRACK NATIONWIDE DELIVERY

We use our own transport for many deliveries which gives us the flexibility to meet your needs.



SUPERIOR BUILD QUALITY

Our products are built to provide years of service. Using the latest materials and fabrication techniques.



KNOWLEDGEABLE SALES TEAM

Our experienced sales team is always available to help you find the best possible solution.



INNOVATIVE WELL DESIGNED PRODUCTS

Our products are designed to meet the ever changing needs of our customers. We continually provide new solutions.



EXCEPTIONAL VALUE FOR MONEY

Through careful design and production processes you can be assured that our products will provide value for money over their working life.



EXTENSIVE PRODUCT RANGES

To give you the best solution we maintain an extensive range of products for complete flexibility.

CONTENTS

CUSHIONS

3

CUSHIONS

Harvest Cushions

05

System Cushions

08

Ultra-Thin Cushion

17



CUSHIONS

4

INTRODUCING
CUSHIONS

The Harvest range of cushions is intended for patients who are sitting for long periods.

Designed to provide increased comfort, they are also suited to those patients who are susceptible to pressure ulcers and can help prevent ulcers developing.

For an extensive range of needs, Harvest Healthcare have an extensive product range, able to meet your requirements. Most cushions are supplied with a waterproof, vapour permeable, stretch cover- protecting the inner foam. The key differences between each type of cushion is the material used for the inner core.

Foam provides a cost-effective solution for low risk needs and can be engineered to create different characteristics. They are normally used for lower risk applications.

Viscoelastic foams have inherently different properties depending on the type and typically provide better pressure distribution. Often used with a composite base made of other foam materials, this makes them more suited to higher risk applications.

Viscoelastic gels will mould to the body shape to distribute pressure and reduce friction and shear effects. Typically used for higher risk applications, they have an added benefit of being impermeable to water and can be sterilised.



t +44 (0)1709 377172 f +44 (0)1709 377173

e sales@harvesthealthcare.co.uk

www.harvesthealthcare.co.uk

MODULAR FOAM

CUSHIONS

5

CUSHIONS

MODULAR FOAM CUSHION

The foam selected is outstanding; a high resilience reflex grade, to stabilise and support, with castellations to distribute pressure and allow air circulation.

It has a navy blue, two-way stretch cover with antislip base.

LOW MED HIGH V HIGH
PRESSURE RELIEF RISK




SIZE OPTIONS

MAX WEIGHT KG

HFC121 - 430mm x 430mm x 95mm	95
HFC116 - 500mm x 500mm x 95mm	108

PRESSURE MAPPING



70kg patient



t +44 (0)1709 377172 f +44 (0)1709 377173
e sales@harvesthealthcare.co.uk
www.harvesthealthcare.co.uk

VISCO SANDWICH VISCO TOP

VISCO SANDWICH CUSHION

The Visco Sandwich Cushion is our two-way turn, viscoelastic foam cushion.

Layers of viscoelastic foam either side of a standard layer of foam, provide the sandwich effect, which offers the user more support and comfort.


It has a navy blue, two-way stretch cover with antislip base.

LOW MED HIGH V HIGH
PRESSURE RELIEF RISK

SIZE OPTIONS

MAX WEIGHT KG

HFC128 - 430mm x 430mm x 95mm	127
HFC128.1 - 510mm x 510mm x 95mm	139



VISCO TOP CUSHION

A 50mm layer of viscoelastic foam tops the standard foam base, transforming the cushion to a high risk product. The top layer of foam provides extra support and comfort to high risk patients.

It has a navy blue, two-way stretch cover with an antislip base.

SIZE OPTIONS

MAX WEIGHT KG

HFC129 - 430mm x 430mm x 76mm	127
HFC123 - 510mm x 510mm x 76mm	139

REFLECT

CUSHIONS

7

CUSHIONS

REFLECT CUSHION

A layer of castellated, viscoelastic foam rests on a high density foam base. The castellations allow the foam to breathe as well as providing more patient comfort. The base gives resilience and stability.

It has a navy blue, two-way stretch cover with an antislip base.

LOW MED HIGH V HIGH
PRESSURE RELIEF RISK




SIZE OPTIONS

MAX WEIGHT KG

HFC130 - 430mm x 430mm x 100mm	127
HFC132 - 500mm x 500mm x 100mm	139

PRESSURE MAPPING



70kg patient



t +44 (0)1709 377172 f +44 (0)1709 377173
e sales@harvesthealthcare.co.uk
www.harvesthealthcare.co.uk

GEL CUSHION

The Gel Cushion is made of high resilience, impermeable, viscoelastic polyurethane gel with highly elastic properties.

It is designed to enhance existing support surfaces.



LOW MED HIGH **V HIGH**
PRESSURE RELIEF RISK

PRODUCT CODE	SIZE	MAX WEIGHT KG
P311C40401HF	400mm x 400mm	120
P311C45431HF	450mm x 430mm	140

MIXTE

CUSHIONS

9

CUSHIONS

MIXTE CUSHION

Combination of gel and foam;

Reticulated polyurethane gel avoids migration and thus prevents penetration (punching) of the cushion while assisting pressure distribution.

The thick high resilience foam base increases the contact area, providing better distribution, comfort and stability.

The two layers are permanently bonded together - eliminating the risk of separation.

ADDITIONAL FEATURES



Combination of Gel and HR Foam




Visco Gel



40Kg/m3 HR Foam



100% Impermeable



LOW MED HIGH V HIGH
PRESSURE RELIEF RISK

PRODUCT CODE	SIZE	MAX WEIGHT KG
P314C40401HF	400mm x 400mm x 60mm	90kg

t +44 (0)1709 377172 f +44 (0)1709 377173
e sales@harvesthealthcare.co.uk
www.harvesthealthcare.co.uk

DUOFORM

ADDITIONAL FEATURES



Highly
Breathable
Cover



100%
Impermeable



40Kg/
m3 HR
Foam



Combination
of Gel and HR
foam



Multi-
bearing
Shape
Concept



High -
Reticulation
Viscoelastic Gel



Fluid
Gel



DUOFORM®

The Duoform® cushion is an anatomically moulded base of high resilience foam surmounted by a two part viscoelastic gel. The white pad fluid gel ensures pressure reduction on high-risk areas, while the blue viscoelastic gel resists fluid migration thus preventing penetration and providing added support.

The high-resilience foam base enables the cushion to better adapt to bony protuberances, spreads contact surfaces and gives better pressure distribution.

It has a highly breathable cover with an antislip base.

LOW MED HIGH V HIGH
PRESSURE RELIEF RISK

DUOFORM® CUSHION

PRODUCT CODE	SIZE	MAX WEIGHT KG
P341C42421HF	420mm x 420mm x 80mm	110
P341C45421HF	450mm x 420mm x 80mm	120
P341C51461HF	510mm x 460mm x 80mm	160

t +44 (0)1709 377172 f +44 (0)1709 377173

e sales@harvesthealthcare.co.uk

www.harvesthealthcare.co.uk

VISCOFLEX PLUS

CUSHIONS

11

CUSHIONS

VISCOFLEX® PLUS CUSHION

The Viscoflex Plus cushion is a very high density viscoelastic, moulded cushion with an ergonomically shaped load bearing core.

The core is shaped to give optimal immersion in sacral/ ischial areas, reducing pressure and facilitating blood circulation. The outer foam is moulded to promote stability and increase the contact area. The skin effect of the moulded foam also gives protection from tearing or crumbling.

This ultra-light cushion is specially designed to provide high levels of stability and pressure relief for the user. It is supplied with a vapour permeable, waterproof cover.

ADDITIONAL FEATURES



Highly
Breathable
Cover



Surface
including
Vents



80Kg/
m3 Visco
Foam



LOW MED HIGH **V HIGH**

PRESSURE RELIEF RISK

VISCOFLEX® PLUS

PRODUCT CODE	SIZE	MAX WEIGHT KG
P361CAF42421HF	420mm x 420mm x 80mm	100
P361CA51461HW	510mm x 460mm x 80mm	150

t +44 (0)1709 377172 f +44 (0)1709 377173
e sales@harvesthealthcare.co.uk
www.harvesthealthcare.co.uk

VISCOFLEX PLUS BARIATRIC



VISCOFLEX® PLUS BARIATRIC CUSHION

The innovative Viscoflex® Plus cushion consists of a moulded Viscoelastic foam with an ergonomic foam insert.


This cushion is specially designed to provide high levels of stability and pressure relief for the user. Supplied with a vapour permeable, waterproof cover with an antislip base.

LOW MED HIGH **V HIGH**
PRESSURE RELIEF RISK


VISCOFLEX® PLUS

PRODUCT CODE	SIZE	MAX WEIGHT KG
P361CAF56461HW	560mm x 460mm x 80mm	230
P361CA61501HW	610mm x 500mm x 80mm	270


WITH STANDARD CUSHION



WITH VISCOFLEX®



**ONE
REDUCE THE LOAD
ON RISK ZONES**
Protected risk zones (sacrum/ischial), even in the cases of pelvic retroversion, in order to provide partial relief of bony protuberances through support.



**TWO
ANTI-PUNCHING
SYSTEM**
Specifically positioned in relation to the ischia. Distributes the pressure to the buttock peripheral zone, stabilizing pelvis support and preventing it from reaching the bottom zone.



**THREE
PRE-ISCHIAL
RIDGE**
Blocks the ischia, thus preventing forward sliding.



t +44 (0)1709 377172 **f** +44 (0)1709 377173
e sales@harvesthealthcare.co.uk
www.harvesthealthcare.co.uk

*Please note information supplied here is a guide. If you would like more information or help with making the correct choice, please call our team today.

VISCOFLEX PLUS WITH POMMEL

CUSHIONS

13

CUSHIONS

VISCOFLEX® PLUS CUSHION WITH POMMEL

The Viscoflex® Plus cushion with pommel consists of a Viscoflex® Plus cushion with an innovative stabilising and adjustable pommel, which is completely separate from the cushion. Also included is a chair strap preventing the cushion sliding forward.

The anatomically shaped pommel provides a true alternative to restraints, preventing sliding and falls.

Supplied with a vapour permeable, waterproof cover with an antislip base.

The ultra-light pommel has an adjustable depth and is available in three sizes.

ADDITIONAL FEATURES



Highly
Breathable
Cover



Surface
including
Vents



80Kg/
m3 Visco
Foam



Anatomically
Shaped
Pommel



Adjustable
Pommel
Depth



VISCOFLEX® PLUS CUSHION WITH POMMEL

PRODUCT CODE	SIZE	MAX WEIGHT KG
P361CP42421HW	420mm x 420mm x 80mm	100
P361CP45421HW	450mm x 420mm x 80mm	110
P361CP45461HW	450mm x 460mm x 80mm	130

LOW MED HIGH **V HIGH**
PRESSURE RELIEF RISK

t +44 (0)1709 377172 f +44 (0)1709 377173
e sales@harvesthealthcare.co.uk
www.harvesthealthcare.co.uk

POLYAIR 60

ADDITIONAL FEATURES



All Polyair® cushions feature dual feature compartments which act independently from each other. This design is useful in rectifying pelvic imbalances while offering greater stability and comfort.




The Polyair 60 is best suited to smaller, somewhat active patients independently able to shift their weight and requiring little or no assistance with transfer



The honeycomb layout optimises the density of the air cells and improves contact surface and pressure distribution.



Simple inflation method, with inflation matrix and pressure gauge manometer.



LOW MED HIGH **V HIGH**
PRESSURE RELIEF RISK

P301C - POLYAIR® 60 CUSHION

PRODUCT CODE	SIZE	MAX WEIGHT KG
P301C40401HW	400mm x 400mm x 60mm	80
P301C40431HW	400mm x 430mm x 60mm	90
P301C43431HW	430mm x 430mm x 60mm	110
P301C51501HW	510mm x 500mm x 60mm	180

POLYAIR® 60 CUSHION

The Polyair® cushion is an adjustable air inflated pressure relief cushion.

The air cells are made with Laxprene which is a high durability compound with excellent tear resistant properties, arranged in a honeycomb design. This greatly reduces pressure through an increased surface area.

The cushion adapts to the patient's body form and cell independence and mobility counteracts shearing compartments.

Comes complete with a vapour permeable, waterproof cover.

POLYAIR 100

CUSHIONS

15

CUSHIONS

POLYAIR® 100 CUSHION

The Polyair 100 is best suited to larger, less active patients with very limited mobility who are at a very high risk of developing, or have already developed, pressure sores.




All Polyair® cushions feature dual feature compartments which act independently from each other. This design is useful in rectifying pelvic imbalances while offering greater stability and comfort.




The honeycomb layout optimises the density of the air cells and improves contact surface and pressure distribution.

CUSHION INFLATED BY HAND



Under Inflated



Over Inflated




Best suited to larger, less active patients with very limited mobility.



Simple inflation method, with inflation matrix and pressure gauge manometer.

INFLATED POLYAIR® CUSHION



Efficient Prevention



LOW MED HIGH V HIGH
PRESSURE RELIEF RISK

P302C - POLYAIR® 100

PRODUCT CODE	SIZE	MAX WEIGHT KG
P302C40401HW	400mm x 400mm x 80mm	80
P302C40461HW	400mm x 460mm x 80mm	120
P302C51461HW	510mm x 460mm x 80mm	160
P302C51501HW	510mm x 500mm x 80mm	180

t +44 (0)1709 377172 f +44 (0)1709 377173
e sales@harvesthealthcare.co.uk
www.harvesthealthcare.co.uk

BACK SUPPORT AND THERAPEUTIC COMBI-PACKS

ADDITIONAL FEATURES



Universal
Attachment
System



80Kg/m3
Visco Foam



BACK SUPPORT & THERAPEUTIC COMBI-PACKS

This back support helps to improve the position and stability of the user. The support has a multi-bearing surface to enable better pressure relief and comfort. It is made of viscoelastic foam with memory effect.

The cover consists of a two way stretch vapour permeable fabric combined with a high abrasion resistant fabric on areas subject to greater body contact.


Available in three sizes.

BACK SUPPORT CAN BE USED TO:


Improve positioning in a wheelchair

Reduce forward sliding.

REDUCTION OF PRESSURE POINTS



Without Back
Support Cushion



With Support
Cushion

P951D- BACK SUPPORT & THERAPEUTIC COMBI-PACKS

PRODUCT CODE	SIZE	WHEELCHAIR WIDTH
P951DS1HW	440mm x 390mm x 90mm	420mm - 500mm
P951D1HW	490mm x 440mm x 90mm	360mm - 440mm
P951DL1HW	490mm x 500mm x 100mm	300mm - 380mm

t +44 (0)1709 377172 f +44 (0)1709 377173

e sales@harvesthealthcare.co.uk

www.harvesthealthcare.co.uk

ULTRA THIN

CUSHIONS

17

CUSHIONS

ULTRA THIN CUSHION

The Ultra-Thin Cushion has been specifically developed for use with a range of chairs.

Available separately, each 35mm thick cushion provides enhanced pressure distribution to the existing chair cushion, without compromising user posture. The thin nature of the cushion also ensures the patient remains safely in the chair as the arms can still be used as a barrier.


The viscoelastic foam offers extra comfort, support and pressure distribution.

This lightweight cushion is also available as a leg and lumbar system and comes with a full strap system to secure it to the chair.


Comes with a waterproof, vapour permeable cover with an antislip base.



REDUCTION OF PRESSURE POINTS



Without
Cushion Fitted



With
Cushion Fitted



LOW MED HIGH V HIGH
PRESSURE RELIEF RISK

ULTRA THIN CUSHION

PRODUCT CODE	TYPE	SIZE	MAX. PATIENT WEIGHT
HFC143.1	Seat Only	430 x 430 x 35mm	24 Stone / 158kg
HFC143LEG	Leg	440 x 430 x 25mm	24 Stone / 158kg
HFC143LUMB	Lumbar	700 x 430 x 25mm	24 Stone / 158kg
HFC143COM.1	3-Piece Set		24 Stone / 158kg

t +44 (0)1709 377172 f +44 (0)1709 377173

e sales@harvesthealthcare.co.uk

www.harvesthealthcare.co.uk

NOTES

NOTES

NOTES

NOTES



Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

T +44 (0)1709 377172

F +44 (0)1709 377173

E sales@harvesthealthcare.co.uk

www.harvesthealthcare.co.uk