

Pump Specifications

	Harvest 1	Harvest 2	Harvest 3
Cycle Time	10 Minutes	10 Minutes	10/15/20/25 Minutes
Ideal For	Community and Care Homes	Community and Care Homes	All Sectors
Pressure Display	None	LED	LED
Air Flow	7 lpm	9 lpm	12 lpm
Dimensions	100 x 130 x 280mm	130 x 140 x 270mm	240 x 120 x 340mm
Low Pressure Alarm	Audible and Visual	Audible and Visual	Audible and Visual
Static Mode	✓	✓	✓
Power Failure Alarm	✗	✓	✓
Lock Out	✗	✓	✓
Seat Inflation	✗	✗	✓
Maxi Firm Mode	✗	✗	✓