



harvesthealthcare® equipped to care

SLING DIMENSIONS, for guidance only. All measurements in cm

B - BACK DEPTH | W - WIDTH | L - LENGTH | LG - LEG | HS - HEAD SUPPORT | BB - BODY BELT

HSL100 Universal Deluxe	B	W	L	LG	HS	BB
Extra Small	52	77	x	53	x	x
Small	58	86	x	58	x	x
Medium	68	96	x	68	x	x
Large	73	116	x	75	x	x
Extra Large	73	129	x	75	x	x

HSL400 In-Chair Hammock	B	W	L	LG	HS	BB
Extra Small	x	81	85	x	x	x
Small	x	86	91	x	x	x
Medium	x	96	105	x	x	x
Large	x	112	113	x	x	x
Extra Large	x	120	120	x	x	x

HSL200 Toilet / Dress	B	W	L	LG	HS	BB
Extra Small	28	60	92	64	x	114
Small	31	70	98	67	x	122
Medium	37	82	119	82	x	132
Large	39	101	126	87	x	142
Extra Large	41	111	130	89	x	152

HSL500 Stand Aid	B	W	L	LG	HS	BB
Extra Small	30	76	x	x	x	*114
Small	33	86	x	x	x	*122
Medium	36	92	x	x	x	*132
Large	40	98	x	x	x	*142
Extra Large	43	104	x	x	x	*152

	Binding	Straps	Loops	Handles
Extra Small	Blue	Blue	Varied	Blue
Small	Yellow	Yellow	Varied	Yellow
Medium	Green	Green	Varied	Green
Large	Red	Red	Varied	Red
Extra Large	Black	Black	Varied	Black
Special	Blue	Blue	**Blue	Blue

* All slings can be identified for size, model, and type of sling by the label.

** Special slings will have blue loops and one coloured loop to identify its size unless it is a complete bespoke product in which case refer to label.



harvesthealthcare®

t +44 (0)1709 377172 **f** +44 (0)1709 377173
e sales@harvesthealthcare.co.uk
www.harvesthealthcare.co.uk