


GRAB HANDLE FOR WOBURN BEDS



Product Code: HLB685
HLB686

Description:

- Compatible with all the Woburn bed range; use HLB686 for the Ultimate bed.
- Assists the patient to sit up and get in and out of bed.
- Fits securely to the bed frame - no additional fixings required.
- Durable epoxy coating.
- Can be used without removing the side rails.

Weight: 3.7 kg

Height: 610mm
(From platform)

Width: 240mm

Safe Working Load: 60kg

Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Grab Handle is guaranteed for one year from date of purchase.

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 2: Nov 2017