

AIR-TECH ADVANCED HEEL PROTECTION HYBRID MATTRESS



Product Code:

HHMAIRT002

Description:

- The Air-Tech Advanced Heel Protection mattress combines the comfort of a static mattress with the pressure relief of alternating cells
- Foam within cells enables the mattress to be stepped up or down with or without a pump providing flexibility for all your nursing requirements.
- 20mm 12 MicroCell tapered heel section is designed to provide zero heel pressure in both static and dynamic modes.
- Firm 'U' shaped foam base for stability, roll prevention and to aid transfers
- Castellated soft and resilient central foam insert with zonal areas designed for support and maximum pressure relief
- Cells powered by the quiet and efficient pump
- Fully welded, antibacterial, waterproof, two way stretch cover
- Durable PVC base

Dimensions:

1990 x 900 x 150mm

Mattress Weight:

14kg

Max. User Weight:

247kg/39 Stone (With Pump)

Pressure Relief:

Very High Risk

Pump Technical Specifications

Name:

Harvest 4

Dimensions:

265 x 175 x 90mm

Cycle time:

12 minutes

Number of Cells:

10 Alternating Plus 12 Alternating MicroCells

Material:

PU material

Cell Cross Section Style:

Single / double cell castellated foam insert with a PU foam surround

Cover:

Welded, two way-stretch, vapour permeable, PU coated cover in plum, zipped to a durable PVC coated base section. Silver-based antibacterial properties in the PU coating are wash proof and environmentally friendly. Cover does not contain latex.

Cleaning Guidance:

Cover should be wiped down with hot water and mild detergent.

If heavily soiled, the cover may be machine washed up to 85°C and tumble dried. For disinfection use sodium hypochlorite solution max 1000ppm. Use no bleach, biological or phenolic-based cleaners. Dry completely before storing.

Warranty

In terms of Harvest Healthcare Guarantee and Warranties the Air-Tech Enhanced Heel Protection mattress is guaranteed for three years from date of purchase.

Exclusions:

1. Removable covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded
3. Components of pumps which are recommended for annual service, are guaranteed for one year.

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact.
(Mattress risk categories are to be used as a guide only)

Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 1: Oct. 2018