

EVE- BATH & SHOWER LIFT



Description:

- The one-of-a-kind, modern EVE shower and bath lift is just the right tool, depending on the mobility of your resident, to be able to safely and comfortably shower or bathe.
- The EVE is a highly manoeuvrable and easy-to-move lift that is versatile and can be used in any existing care bathtub.
- The battery-powered electric motor makes for an infinitely variable and quiet height adjustment that will not agitate the resident. With the associated handheld remote control, the lift can be operated simply and safely by one staff member.
- An integrated service module reports technical problems with the lift as well as necessary maintenance.

Safe Working Load:

185kg (29 Stone)

Lifting Range:

495mm - 1065mm

Key Features:

- Made entirely of aluminium and stainless steel
- 100% functional as a shower chair
- Simple operation by only one staff member
- Fits any modern care bathtub
- Unrestricted transfer from the Ascent Pro to the seat lifter
- The carrier frame is formed in a way that allows the resident to sit down on the seat without hindering the legs and heels
- Ergonomically designed safety steering handles ensure accurate manoeuvring
- Low entry height of only 49.5 cm
- Closed, easy-to-manoeuve double castors
- Handheld control unit with visual indication of battery level
- Battery and charger are compatible with all Harvest mobile lifts



Double Castors
Closed, easy-to-manoeuvre double castors, two with brakes.



Optional Scale
Optional: Integrated scale powered by the main battery.



Comfort Seat
Comfortably shaped seat with hygiene opening.



Fits any Care Bathtub
The EVE fits in any modern care bathtub.



Aluminium Lifting Column and Carrier Frame
Aluminium components make the lift extremely lightweight.



Safety Steering Handles
Enclose the resident and ensure accurate manoeuvring.

Technical Specifications

Lifting Power:	570mm
Length:	950mm
Width:	680mm
Height:	1,240mm
Entry Height (min):	495mm
2 Double Castors (Front):	100mm diameter
2 Double Castors (rear, with fixing):	100mm diameter
Free Space under Carrier Frame:	30mm
Carrier Frame Height:	140mm
Battery:	24V

Accessories/ Options

Digital Scale
Replacement Battery
External Wall Charging Station

Supplier Contact Details:

t: +44(0)1709 377172 e: sales@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products.

Products and data continually undergo development and may therefore deviate from the details given.

Version 1: May 2019