



MARINA SHOWER TROLLEY

**Description:**

- The Harvest Healthcare Marina Shower Trolley allows for effective hygiene and cleaning of immobile and bed bound patients.
- Foam padded PVC liner with side supports and end panels for comfort and security
- Supplied with foam padded PVC patient support cushion
- Fold down side support rails to assist patient transfer

Adjustable Height:

410mm-940mm

Safe Working Load:

37 Stone / 240kg

Overall Size:

1950mm x 740mm

Features:

Low gradient tilt towards foot end for water drainage via integral waste outlet
Plug and drain hose supplied complete with storage bracket
Bumpers fitted to all corners for protection during transport procedures
Removable bumpers to reduce the width to 65cm to access tight doorways
Lowers to wheelchair height for ease of patient transfer
Purpose made precision bearings on all pivot points
Clearance beneath frame for mobile hoist

Supplier Contact Details:

t: +44(0)1709 377172 e: sales@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products.
Products and data continually undergo development and may therefore deviate from the details given.

Version 1: Jan 2019