

QUICK MOVE



Product Code: ROMS400801334

- Description:**
- Developed specifically for users with impaired balance, the QuickMove allows users to be transferred with the aid of a carer.
 - QuickMove provides a more secure sit-to-stand procedure for both user and carer, and, at the same time, it contributes to safeguarding and strengthening the user's abilities.
 - Suitable both for long-term and for occasional needs.
 - Made of non-corrosive materials.
 - Compact and convenient and very easy to manoeuvre which makes it very useful in most spaces and situations.

Safe Working Load: 170kg (27 Stone)

Length: 980mm

Width: 605mm

Height: 1120mm

Base Width: 510-900mm

Features: Adjustable base width which makes it easier to get close and in the right position for the user
 Adjustable buttock cushions can be angled to make room for the user and caregiver
 Adjustable lower leg support
 Support of the cushions behind the buttocks and the knee-brace facilitates a very stable and safe transfer

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 1: Jan 2019