



ASCENT PRO 200 STAND AID



Description:

- The Ascent Pro 200 Stand aid is ideal for patients with reduced core or leg strength to give them the extra assistance they need to enable them to reach a standing position in an efficient manner.
- It has been designed with client comfort in mind, with easy to grip handles adjustable to suit a range of heights, and soft padded leg supports with integrated shock absorbers.
- The Ascent Pro is battery operated for full manoeuvrability.
- The overload protection and emergency lowering systems ensure that a patient will always be gently lowered to a seating position in the event of power failure.
- It also boasts battery operated spreadable legs to offer a stable lifting position without impacting on manoeuvrability.

Safe Working Load: 200kg (31.5 stone)

Adjustable Height: 1168mm-1781mm

Length: 1027mm

Width: 700/1044 (closed/spread)

Key Features: Battery operated height and leg adjustment (24 volt battery)
Handset operation
Manual emergency lowering system
External wall charging station
Flexible extra soft knee supports with safety straps
4 Castors (2 double and 2 lockable)
Ergonomic handles
5 Year warranty

Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Ascent Pro Stand Aid is guaranteed for five years from date of purchase.

Supplier Contact Details:

t: +44(0)1709 377172 e: sales@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products.
Products and data continually undergo development and may therefore deviate from the details given.

Version 1: Mar 2019