

# ELITA LOW PROFILING BED



**Product Code:**

HLB559

**Description:**

- The Elita Low is a premium profiling bed which lowers to just 230mm from the floor. Ideal for patients at risk of bed falls the Elita Low provides exceptional care whilst maintaining a warm aesthetic.
- Four profiling sections allow for optimum comfort and clinical support, and can be independently adjusted by users using the ergonomic lockable handset.
- Fully customisable, the Elita Low comes with multiple colour options, wraparound head and foot boards, and optional Trendelenburg and auto regression.

**Safe Working Load:**

175kg / 27 Stone

**Max. Patient Weight:**

135kg / 21 Stone

**Dimensions:**

1020 x 2150mm

**Adjustable Height:**

230mm - 630mm

**Angle of Back Rest:**

70 Degrees

**Angle of Thigh Rest:**

20 Degrees

**Key Information:**

Lockable handset  
 Electric backrest and knee break with mechanical lower leg support  
 Auto regression (optional)  
 Trendelenburg (optional)  
 Stylish wooden side rails Wraparound head and foot boards  
 Available in 8 colours  
 4 braked castors  
 Comes with a 5-year warranty

**Supplier Contact Details:**

t: +44(0)1709 377 172

e: enquiries@harvesthealthcare.co.uk

w: www.harvesthealthcare.co.uk

Sheaf House, Bradmarsh Way, Bradmarsh Business Park, Rotherham, S60 1BW

## Available Colours:



LIGHT BEECH R 5311



MAPLE R 5703



BEECH R 5111



ACACIA R 5052



CHERRY R 4991



LIGHT OAK R 4272



NABUCCO R 3144 L



WALNUT H 3704

Wood swatch samples available upon request.

## Warranty

In terms of Harvest Healthcare Guarantee and Warranties the Elita Low Profiling Bed is guaranteed for five years from date of purchase.

### Exclusions

1. Actuators are guaranteed for 2 years
2. Handsets are guaranteed for 6 months

### Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products.

Products and data continually undergo development and they may therefore deviate from the details given.

**Version 1: June 2018**