

VISCOFLEX PLUS CUSHION



Product Code:	P361CAF42421HF
Description:	<ul style="list-style-type: none"> • A very high density viscoelastic, moulded cushion with ergonomically shaped load bearing core. • Core is shaped to give optimal immersion in sacral/ ischial areas. • Outer foam moulded to promote stability and increase contact area. • Available in a range of size and weight capacities.
Dimensions:	420 x 420 x 80mm
Weight of Cushion:	1.25 kg
User Weight:	40-100 kg
Pressure Relief Risk:	High to Very High
Foam Density:	
Core:	40 kg/m ³
Outer	80 kg/m ³
Cover:	Anti-allergenic, antimicrobial, antibacterial, anti-fungal, two way stretch, vapour permeable, washable, welded seams, two side zip. Printed user instructions on cover.
Fire Retardancy:	NF EN ISO 597 - 1&2

Cleaning Guidance:

The cushion foam may be hand washed with mild detergent. Do not immerse in water. Allow to dry naturally. Cover should be wiped down with hot water and mild detergent or machine washed up to 95°C and tumble dried. For a comprehensive list of tested cleaning products see our website.

Warranty

In terms of Harvest Healthcare Guarantee and Warranties the Viscoflex Plus cushion is guaranteed for three years from date of purchase. Exclusions:

Exclusions:

1. Removable cushion covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

**Important Safety Notice**

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact.
(Mattress risk categories are to be used as a guide only)

Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 2: Oct 2017