



# WOBURN STANDARD BED & ACCESSORIES




**HLDU582/HLDU582.14 - ULTIMATE LIFTING POLE** COMPATIBLE WITH OUR EXCLUSIVE WOBRUN ULTIMATE BED, THIS ULTIMATE LIFTING POLE ASSISTS THE USER TO SIT UPRIGHT IN BED.



**HLDW300/HLDU300 - BED HOOKS** FOR USE WITH ANY PUMPS, WHICH POWER AN ACTIVE MATTRESS THE HOOKS ATTACH TO THE BED END, AND CAN THEN HOLD THE PUMPS.



**HJR075 - BED ALARM SYSTEM CONTROL BOX** IS DESIGNED FOR USE WITHIN A BED, NOTIFYING CARERS OF PATIENT MOVEMENT. THIS HELPS REDUCE THE LIKELIHOOD OF FALLS.



**HLDW / HLDWC / HLD - BED LENGTH EXTENSION KIT** IS USED TO EXTEND YOUR PROFILING BED BY AN ADDITIONAL 200MM. SOME USERS ARE NOT COMFORTABLE IN A TYPICAL 2000MM LONG PROFILING BED, WHICH CAN CAUSE ENTRAPMENT ISSUES