

ASCENT EXTRA 225 STAND AID



Product Code: HLA225

Description:

- The Ascent Extra 225 Stand Aid has been developed in conjunction with healthcare professionals.
- This stand aid features excellent build quality, functionality and attractive looks with ease of use and simplicity.

Safe Working Load: 225kg (35 Stone)

Adjustable Height: 887mm - 1740mm

Length: 1215mm

Width: 1330mm

Key Features:

- Mounted on 4 lightweight castors allowing for 360 swivelling for turning and rolling movement, the rear 2 castors are braked.
- Made of steel with a silver powder coat
- Linak Jumbo Care control system with LCD display
- Designed to reduce infection control issues
- High speed, soft stop/start Linak actuator reduces transfer time.
- Removable non slip foot plate.
- Height and angle adjustable knee pad ensures patient comfort.
- Electric leg adjustment and a neoprene calf strap to prevent tissue damage
- CE marked
- Optional transfer seat available

Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Lifter is guaranteed for five years from date of purchase.

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products.
Products and data continually undergo development and may therefore deviate from the details given.

Version 1: Dec 2018