



harvesthealthcare®


equipped to care

HLDW100 UNDERBED LIGHT WITH SENSOR

The Underbed Light is fitted with a sensor that will detect movement when in close proximity to the bed. This provides the user with a clear light that illuminates the floor around the bed, making it easier to navigate in a dark room or at night.

This kit provides all that you need to install the Underbed Light. The light is positioned directly onto the bed frame and connects to the bed's own power supply.

Kit Contents



1- Mounting Clips

2- Light Tube

3- Sensor

4- Y Connector

5- Cable Clips



Sheaf House, Bradmarsh Business Park, Bradmarsh Way, Rotherham S60 1BW
T +44 (0)1709 377172 F +44 (0)1709 377173 E sales@harvesthealthcare.co.uk

