

3-LAYER VISCO TOP CUSHION



Product Code:

HFC129.01

Description:

- A high risk viscoelastic and foam bonded cushion with improved immersion of visco top and a combination of a second layer of viscoelastic foam and a firm foam base add support.
- The foam combination enhances pressure distribution and combines stability and comfort for the patient at high risk.
- Two way stretch cover with antislip base.
- Cover is latex free and clearly printed, including cleaning instructions.

Dimensions:

430 x 430 x 75mm

Max. User Weight:

158kg/25 stone

Pressure Relief Risk:

High Risk

**Cover:
(Available Separately)**


Four-way stretch, vapour permeable, PU coated top in navy blue. Zip fastened on 2 sides with single flap cover to prevent fluid ingress. Base section is antislip PVC to reduce unwanted movement.

**Foam:
Density and Hardness**

Top: 58-62kg/m³ 70-100N
Middle: 58-62kg/m³ 100-130N
Base: 34-37kg/m³, 90-120N

Fire Retardancy:

BS 5852-1: 1988 crib 5



Visco Top Foam: 58-62kg/m³ 70-100N

Visco Foam: 58-62kg/m³ 100-130N

Firm Foam Base: 34-37kg/m³, 90-120N

Cleaning Guidance:

Wipe down cover with hot water and mild detergent. Or machine wash and dry. Should the foam interior become damp remove cover and allow to dry naturally. Use no bleach. Dry completely before storage. For a comprehensive list of tested cleaning products see our website.

Warranty

In terms of Harvest Healthcare Guarantee and Warranties the 3-Layer Visco Top Cushion is guaranteed for three years from date of purchase with the following exclusions:

Exclusions:

1. Removable cushion covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact.
(Mattress risk categories are to be used as a guide only)

Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 1: Oct 2018