

OLYMPIC PATIENT TURNER



Code: HLA045

Description:

- Harvest Healthcare's Olympic Turner is designed with the comfort, safety and ergonomics of both caregiver and user in mind
- Four bar handle makes it suitable for users of different heights
- Contoured, padded dual knee pads are height and width adjustable for comfort
- Four castors and two central wheels to enable effortless transfers over short distances
- Standing platform covered with a non-slip mat for safety and is low in height to fit under most divan beds to make it easy to step onto
- Suitable for use with those who can support their own weight
- Users pull themselves into position using their own strength and muscles, helping them to stay active
- Optional safety belt available in 3 sizes
- Easy to assemble, no tools required
- Plastic protective bumpers
- Cutout fits around toilet

Safe Working Load: 160kg (25 Stone)

Front Castors: 35mm

Rear Castors 125mm

Assembled Unit: 19kg

Key Features:

- Helps to assist a user with sit to stand transfers
- Lightweight and easy to manoeuvre
- Padded leg support for added comfort
- Multiple grip options for both user and carer
- Olympic Strap available to enable a carer to assist a user to their feet
- Independent locking rear brakes
- 3600 turning radius

Optional Extras:



Freeway Safety Belt

- Provides additional lumbar support for the user
- VELCRO® Brand straps and the adjustable buckle fastening support belt and provide additional security
- Can be used by one or two carers to assist the user to stand
- Available in S, M, L
- Product Codes: 2002ACCRAISBELTS, 2002ACCRAISBELTM and, 2002ACCRAISBELTL



Freeway Ladder Belt

- Sliding sleeve provides comfort for the user and ease of use for the carer
- Can be used by one or two carers to assist a user to their feet
- Dimensions: 590(W) x 120(D)mm
- Product Code: 2002ACCRAISLAD

Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Olympic Turner is guaranteed for one year from date of purchase.

Supplier Contact Details:

t: +44(0)1709 377172 e: sales@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products.

Products and data continually undergo development and may therefore deviate from the details given.

MAY 2021