



Moving & Handling Product Range

POWER TO DELIVER

WHO WE ARE & WHAT WE DO

Harvest Healthcare has more than 20 years' experience in design and manufacture of pressure care equipment, beds, moving, handling and showering aids.

We offer a variety of support services from mattress repair and decontamination, service and repair provision and client assessment through to project management for installed solutions.

SUPERIOR PRODUCT SOLUTIONS

We provide a complete pressure care equipment, beds and moving & handling product portfolio consisting of static and active (dynamic) mattresses and cushions, profiling beds, hoists, ceiling track hoists, mobile hoists, stand aids, gantry systems, slings, shower chairs, shower stretchers and much more.

HARVEST HEALTHCARE YOUR FIRST CHOICE...



TOTAL SUPPORT TEAM

We offer a full package to support you from installation to service and maintenance.



HARVEST FAST TRACK NATIONWIDE DELIVERY

We use our own transport for many deliveries which gives us the flexibility to meet your needs.



SUPERIOR BUILD QUALITY

Our products are built to provide years of service by using modern materials and fabrication techniques.



KNOWLEDGEABLE SALES TEAM

Our experienced sales team is always available to help you find the best possible solution.



INNOVATIVE WELL DESIGNED PRODUCTS

Our products are designed to meet the ever changing needs of our customers. We continually provide new solutions.



EXCEPTIONAL VALUE FOR MONEY

Through careful design and production processes you can be assured that our products will provide value for money throughout their working life.



EXTENSIVE PRODUCT RANGES

To give you the best solution we maintain an extensive range of products for complete flexibility.

CONTENTS

1. What Our Service Offers You	04
2. Your Contacts	05
3. Stand Aids, Hoists and Transfer Aids	06
4. Slings	18
5. Bed Management Systems	24
6. Showering and Toilet Aids	28

WHAT OUR SERVICE OFFERS YOU

- A tailored service to suit your requirements.
- A skilled, experienced and professional workforce.
- A cost effective solution.
- Fully traceable work program and database with comprehensive reporting.
- A personal service with named contacts who will look after your needs.
- A complete parts service.



PLANNED/PREVENTATIVE MAINTENANCE

- Lifting Equipment – examination and servicing (LOLER).
- Beds – Inspection and service (PUWER).
- Mattresses – Inspection, decontamination and service/repairs.
- Wheelchairs – Inspection and service (PUWER).
- Suction machines/Nebulisers – Inspection and Service.
- Training packages for your maintenance staff.
- Production and maintenance of your asset register.

INSTALLATION & TRAINING

We have a team of highly experienced installation engineers who are available to install and demonstrate our products for you. We also provide moving and handling advice and training for you to be confident in using our products both efficiently and safely.

MAINTENANCE & REPAIRS

Our focus is to provide you with a first-time, quick and efficient fix, whenever possible. As a manufacturer of healthcare equipment we have a comprehensively stocked parts department.

REACTIVE MAINTENANCE (ALL PRODUCTS)

- On or off-site repair options.
- Same day 24/7 response time.
- Mattress management – repairs, decontamination, exchange.
- Training packages for HH products.
- Telephone support for troubleshooting and spares.

Whether you want us to support you when you require, or to take on a complete asset management role, give us a call and let us know how you would like us to tailor a package for you.

MATTRESS DECONTAMINATION

Our air mattress repairs and decontamination division operates an experienced and dedicated service in repair, maintenance and cleansing of most makes of active pressure relief systems. Our on-site laundry uses OTEX technology which injects OZONE into the wash cycle destroying all bacteria assuring that your equipment is certified decontaminated with full traceability.

EQUIPMENT MANAGEMENT

From hospitals to walk-in centres and rehabilitation clinics, we have a dedicated department that specialises in the maintenance and servicing of equipment.

HERE TO HELP

Harvest Healthcare Contact Number: 01709 377172


Option 1 - PRODUCT SUPPLY – sales@harvesthealthcare.co.uk

Option 2 - SERVICE/REPAIR SCHEDULING – servicing@harvesthealthcare.co.uk

Option 3 - MATTRESS REPAIR DEPARTMENT – mattress.repairs@harvesthealthcare.co.uk

Option 4 - ACCOUNTS – accounts@harvesthealthcare.co.uk

FOR MORE INFORMATION VISIT **WWW.HARVESTHEALTHCARE.CO.UK**



FIND US TODAY

Harvest Healthcare Ltd
Sheaf House, Bradmarsh Business
Park,
Bradmarsh Way,
Rotherham S60 1BW

T +44 (0)1709 377172

F +44 (0)1709 377173

E sales@harvesthealthcare.co.uk

www.harvesthealthcare.co.uk

OLYMPIC PATIENT TURNER



harvesthealthcare®
a prismhealthcare company

Harvest Healthcare’s Olympic Turner is designed with the comfort, safety and ergonomics of both caregiver and patient in mind.

Designed for patients that have some capacity to pull themselves up but need a helping hand, the Olympic Turner requires patients to use their own strength and muscles, helping them to stay active.

KEY FEATURES & BENEFITS

- Helps to assist a user with sit to stand transfers
- Lightweight and easy to manoeuvre
- Padded leg support for added comfort
- Multiple grip options for both user and carer
- Freeway Raiser Strap available to enable a carer to assist a user to their feet
- Compatible with the Freeway safety and ladder belts
- Independent locking rear brakes
- 360° turning radius

1

year warranty

160kg/25st
max user weight



2002ACCRAISLAD	Freeway Ladder Belt
2002ACCRAISBELTS	Freeway Safety Belt - Small
2002ACCRAISBELTM	Freeway Safety Belt - Medium
2002ACCRAISBELTL	Freeway Safety Belt - Large



Freeway Safety Belt



Freeway Ladder Belt
590(W) x 120(D)mm

Technical Information

Product Code	Overall Height	Assembled Unit Weight	Front and Rear Castors/Wheels
HLA045	1201mm	19kg	35mm/125mm

OLYMPIC PATIENT MOVER



harvesthealthcare®
a prismhealthcare company

Simple and intuitive to use, the Olympic Patient Mover is an attractive universal resident transport solution. Using the Olympic Patient Mover allows resident transfers with minimal handling, supports residents to move comfortably from a seated to a standing position, and promotes early mobilisation.

KEY FEATURES & BENEFITS

- Provides easy transportation to the bathroom and toilet
- Offers a suitable alternative to a wheelchair
- Allows for better eye contact between the user and carer
- Padded knee support for extra comfort
- Manual leg open & closing
- Compatible with the Mackworth back support belt


1

year warranty

180kg/28st
max user weight



1201TAB510C	Mackworth Back Support Belt S
1201TAB520C	Mackworth Back Support Belt M
1201TAB530C	Mackworth Back Support Belt L
1201TAB540C	Mackworth Back Support Belt XL



Mackworth
Back Support Belt

Technical Information

Product Code	Length	Width	External Width (Legs Open)	Turning Radius
HLA046	1500mm	650mm	1060mm	1125mm

QUICKMOVE 170 ADJUSTABLE PATIENT MOVER



The QuickMove supports the user when rising, using their own strength, to a standing position and during transfers.

The QuickMove has been developed for patients with impaired balance and strength in their legs. The user must have some strength in the arms and must be able to bear weight to actively participate in the transfer.

KEY FEATURES & BENEFITS

- Small and convenient
- Stimulates natural movement pattern
- Adjustable seat cushions that can be angled
- Adjustable lower leg support/knee-brace
- Adjustable base width

3 year warranty



170kg/27st
max user weight



Technical Information

Product Code	Length	Width	Height
ROMS400801334	980mm	605mm	1120mm

ASCENT 180 STAND AID



The Ascent 180 Stand Aid has an impressive range of ergonomic features to provide a safe and dignified lift.

It is uniquely designed to provide three functions; as a transfer hoist, as a sit to stand hoist and as an aid to rehabilitation by assisted standing.

KEY FEATURES & BENEFITS

- Ability to lift a wide range of people from 4ft – 6ft 5"
- Compact and easily manoeuvred
- Profiled push handles designed in conjunction with the REBA assessment tool allows for safe moving and positioning
- Powered leg opening allows for a stable and safer lift without restricting access around chairs and beds
- Front castors and rear castors for maximum control and positioning on all floor surfaces
- Multi-directional knee pad allows for individual adjustment to provide full support
- Rear foot bar helps to negotiate the stand aid over uneven floor surfaces and door thresholds

5

year warranty

180kg/28st
max user weight



For Clip Stand Aid Slings
HLA038-Manual Leg
HLA039-Powered Leg



For Loop Stand Aid Slings
HLA041-Manual Leg
HLA042-Powered Leg

Technical Information

Product Code	Length	Height (Lowest/Highest)	Leg Inside Width (Closed/Spread)	Lifting Range
HLA038-Clip with Manual Leg HLA039-Clip with Powered Leg HLA041-Loop with Manual Leg HLA042-Loop with Powered Leg	119mm	1145mm/1720mm	540mm/1010mm	885mm

ASCENT PRO 200 PREMIUM STAND AID



With a sleek, ergonomic design the Ascent Pro 200 provides ideal assistance for transfers while maintaining a home from home aesthetic.

Suitable for transferring patients into a wheelchair, onto toileting facilities, onto a bed, or for many other assisted transfers the Ascent Pro 200 is also excellent for standing exercises and promoting mobility.

KEY FEATURES & BENEFITS

- Premium, sleek, aesthetic design
- Battery operated height and leg adjustment
- Soft stop/start motion
- Handset operation with battery level display
- LCD display for service notifications, actuator cycles and overloads
- Ergonomic handles
- Manual emergency lowering system
- Flexible extra soft knee supports with safety straps
- Removable non-slip foot plate
- Enclosed castors keep out dirt and are easy to maintain

5 year warranty



200kg/31.5st
max user weight



OPTIONAL EXTRAS:

- Spare battery
- Wall charging station

Technical Information

Product Code	Length	External Width (Closed/Spread)	Adjustable Height
HLA006	1027mm	700mm/1044mm	1168mm-1781mm

ASCENT EXTRA 225 PLUS SIZE STAND AID



harvesthealthcare®
a prismhealthcare company

The Ascent Extra 225 stand aid is designed to make everyday mobility tasks simple, safe and stress free.

The ergonomic design allows operation by a single caregiver, and the high speed, soft stop/start motion makes transfers smooth and reduces transfer time. Suitable for transferring residents into a wheelchair, onto toileting facilities, on to a bed or for many other assisted transfers the Ascent Extra 225 is also excellent for standing exercises and promoting mobility.

KEY FEATURES & BENEFITS

- Battery operated height and leg adjustment
- Handset operation with battery level display
- Manual emergency lowering system
- Flexible extra soft knee supports with safety straps
- Electric leg adjustment
- Neoprene calf strap
- Removable non-slip foot plate
- External wall charging station

5 year warranty



225kg/35st
max user weight



Transfer Seat
HLA225.ST

Technical Information

Product Code	Length	External Width (Legs Open)	Adjustable Height
HLA225	1215mm	1330mm	887mm-1740mm

OLYMPIC 180 HOIST



harvesthealthcare®
a prismhealthcare company

The Olympic 180 Hoist benefits from the rigorous ergonomic design procedures and quality checks. It has been constructed from high quality powder coated steel that will provide years of reliable service. The hoist boasts an impressive 180kg lifting capacity, powered leg opening and market leading Linak® actuators and motors.

KEY FEATURES & BENEFITS

- Cost effective solution without compromising on quality
- Robust powder coated steel frame
- Linak control system
- Electric leg opening as standard
- Red emergency lowering screw
- Easy to position with 4 large swivel castors and ergonomic handle
- Tool free assembly
- Single person installation
- Can lift a patient from the floor
- 75mm front castors and 100mm rear castors
- Available with scales unit incorporating QRS which is a dignified, accurate and convenient method of determining body mass by weighing an individual whilst situated in a sling and suspended by the Olympic 180 hoist

5 year warranty



180kg/28st
max user weight



QRS

Technical Information

Product Code	Length	Height (Highest/Lowest)	Leg Inside Width (Closed/Open)	Lifting Range
HLA034 HLA034QRS with QRS	1265	1335mm/2010mm	490mm/1075mm	1755mm

OLYMPIC PRO 185 PREMIUM HOIST



The Olympic Pro 185 hoist provides a beautifully designed, premium solution for labour intensive resident transfers, allowing powered, single carer use. The powered dynamic positioning system (PDPS) provides precise and stable positioning without the need for manual intervention.

The height, leg spreader and carrier frame are all electronically controlled, and the anti-rotation mechanism prevents unintentional rotational movement during lifting and transferring. The carrier frame can be lowered to the floor to allow picking up residents from the ground. Suited to a wide range of care environments.

KEY FEATURES & BENEFITS

- Aesthetic design, extra light weight
- Emergency stop and emergency lowering
- Carrier frame and lifting column made of aluminium and steel
- Enclosed lockable double castors
- Electrically adjustable height, cradle and spreadable carrier frame
- Electrical Positioning System (EPS)
- Anti-rotation mechanism
- Battery compatible with Ascent Pro 200 stand aid
- Slow decline feature, in case of overloading
- LOLER test reminder light

5 year warranty



185kg/29st
max user weight



OPTIONAL EXTRAS:

- Electronic weighing scale (HLA010)
- Spare battery
- Wall charging station
- Two point yoke frame for loop slings

Technical Information

Product Code	Chassis Length/Height	Chassis Width (Closed/Spread)	Maximum Height	Lifting Range
HLA002	1290mm/125mm	700mm/1180mm	2170mm	1400mm

OLYMPIC EXTRA 230 PLUS SIZE HOIST



The Olympic Extra 230 is a sling hoist specially made for lifting, moving and transferring plus size patients. Due to the special construction and the steel components the Olympic Extra is extra light and easy to move, but still offers the required stability for the care of plus size patients.

We provide compatible slings for our hoists.

KEY FEATURES & BENEFITS

- Aesthetic design
- Universal use, also appropriate for bathing
- Extra light and easy to manoeuvre
- Electrical positioning system (EPS) places the patient exactly without additional force
- Electrically adjustable height, cradle and spreadable carrier frame, for easier patient transfer
- Lifts from the floor
- No rotational movements during lifting and transferring due to anti-rotation mechanism
- Hand control with battery charge status displayed
- In-built electronic scales (optional extra)
- Wall charging station (optional extra)

5 year warranty



230kg/36st
max user weight



Pivot frame
To be used with clip slings



2 point yoke frame
To be used with loop slings

Technical Information

Product Code	Length	External Width (Legs Open)	Min Height / Max Height
HLA007	1400mm	1240mm	1340mm / 2170mm

FREEWAY SCALES UNIT INCORPORATING QRS

FREEWAY

A dignified, accurate and convenient method of determining body mass by weighing an individual whilst situated in a sling and suspended by a hoist.

This product is compatible with the Olympic 180.


KEY FEATURES & BENEFITS

- Fits in-line between the hoist and carry bar
- Load of 10 to 350kg (functional maximum load is limited by the hoist or sling, i.e. 180kg when used with the Olympic 180 hoist).
- QRS (Quick Release System)
- 1 button operation
- Auto shut down after 60 seconds (programmable)
- 15mm LED screen
- Display resolution to 0.2kg
- Class III accuracy

5 year warranty



180kg/28st
max user weight



Technical Information

Product Code	Description
610000	QRS & Scales Assembly
1102CBQRS3HK	3 Hooked Clipped Carry Bar

FREEWAY TRANSACTIVE XTRA CEILING TRACK HOIST




The Freeway Transactive Xtra is a product of world class design and state of the art features, providing safe, reliable and effective hoisting for carers and users alike.

The Transactive Xtra benefits from a low profile design which is unobtrusive, yet provides a much greater lifting height than any other hoist option. Together with the extensive range of track configurations, the Transactive Xtra can provide the perfect solution for any hoisting requirements.

KEY FEATURES & BENEFITS

- Ergonomically designed handset
- State-of-the-art display & diagnostics
- Freeway Spreader Bar as standard
- Emergency stop
- Powered emergency lowering
- Manual lowering & raising facility

2 year warranty



Installed Ceiling Track Systems – multiple options subject to site survey

- Straight Track – Point to Point transfer
- Room Covering – Pick up and transfer to multiple points within the room
- Room to Room – Transfers between multiple linked rooms

Technical Information

Product Code	Description
1001TX270D	Transactive Xtra Drive - 270kg
1001TX270	Transactive Xtra Manual - 270kg
1001TX200D	Transactive Xtra Drive - 200kg
1001TX200	Transactive Xtra Manual - 200kg
1001TX160D	Transactive Xtra Drive - 160kg
1001TX160	Freeway Transactive Xtra Manual - 160kg
1001TX130D	Freeway Transactive Xtra Drive - 130kg
1001TX130	Freeway Transactive Xtra Manual - 130kg

UNIVERSAL DELUXE SLING



This multi-purpose sling is designed to offer excellent support, and to be used for the majority of patients, including for bathing when using super soft material.

Available in polyester which is versatile, breathable and durable.

It is also available in super soft 3-dimensional fabric which is breathable, very durable and provides excellent pressure distribution. Super soft fabric is also fast drying and therefore suitable for bathing.

KEY FEATURES & BENEFITS

- Can be applied and removed while seated in a chair
- Full padded support of upper thigh reduces hip rotation and abduction
- Available in a range of sizes
- Available in Polyester and Super Soft material
- Pocket in the back of the sling to aid proper fitting
- Available with head support
- Also available with clip fastening
- Extra stitching to improve life and safe working load of the sling
- Information stitched into label to protect when washing




300kg/47st (super soft)
max user weight



400kg/62st (Polyester)
max user weight

1 year warranty



Available in
Polyester Fabric



Available in
Super Soft Fabric



Universal Deluxe
Loop Sling
HSL100



Universal Deluxe
Clip Sling
HSL110

Technical Information

Product Code	Sizes	Width Range (cm)	Back Length Range (cm)
HSL100/HSL110	X Small	45-52	48-56
	Small	52-60	53-63
	Medium	60-76	63-73
	Large	76-92	68-78
	X Large	92-100	68-78

IN-CHAIR HAMMOCK SLING



harvesthealthcare®
a prismhealthcare company

The In Chair Hammock Sling is designed to provide full support to the back, neck and shoulders.

Available in polyester which is versatile, breathable and durable.

It is also available in super soft 3-dimensional fabric which is breathable, very durable and provides excellent pressure distribution. Super soft fabric is also fast drying and therefore suitable for bathing.

KEY FEATURES & BENEFITS

- Hip Straps
- Pommel
- Pockets to tuck lifting straps in
- Available in a range of sizes
- Available in Polyester and Super Soft materials
- Also available with clip fastening
- Extra stitching to improve life and safe working load of the sling
- Information stitched into label to protect when washing




300kg/47st (super soft)
max user weight



400kg/62st (Polyester)
max user weight


1 year warranty



Available in
Polyester Fabric



Available in
Super Soft Fabric



In Chair Hammock
Loop Sling
HSL400



In Chair Hammock
Clip Sling
HSL410

Technical Information

Product Code	Sizes	Width Range (cm)	Back Length Range (cm)
HSL400/HSL410	X Small	45-52	48-56
	Small	52-60	53-63
	Medium	60-76	57-67
	Large	76-92	62-72
	X Large	92-100	65-75

TOILET / DRESS SLING



harvesthealthcare®
a prismhealthcare company

The Toilet/Dress Sling has been designed to help facilitate the toileting process, with good access to allow for the removal of clothing.

Available in polyester which makes it versatile, breathable and durable.


KEY FEATURES & BENEFITS

- Easy to fit
- Fully padded
- Adjustable Velcro body belt
- Wide legs for additional comfort
- Slide sheet material on legs to reduce friction and ease fitting
- Available in a range of sizes
- Also available with clip fastening
- Extra stitching to improve life and safe working load of the sling
- Information stitched into label to protect when washing

1

year warranty

400kg/62st
max user weight



Available in Polyester Fabric



Toilet/Dress Sling
HSL200

Technical Information

Product Code	Sizes	Width Range (cm)	Back Length Range (cm)
HSL200	X Small	40-50	26-30
	Small	50-60	29-34
	Medium	62-68	34-40
	Large	68-78	36-42
	X Large	78-85	38-44

STAND AID SLING



The Stand Aid Sling is designed to assist with the standing process using a hoist. It can be used with Harvest Healthcare Ascent Extra Stand Aid.

Available in polyester which makes it versatile, breathable and durable.


KEY FEATURES & BENEFITS

- Easy to fit
- Fully padded
- Adjustable Velcro body belt
- Wide legs for additional comfort
- Slide sheet material on legs to reduce friction and ease fitting
- Available in a range of sizes
- Also available with clip fastening
- Extra stitching to improve life and safe working load of the sling
- Information stitched into label to protect when washing

1

year warranty

400kg/62st
max user weight



Available in
Polyester Fabric



Stand Aid Loop
Sling
HSL500



Stand Aid Clip
Sling
HSL510

Technical Information

Product Code	Sizes	Width Range (cm)	Back Length Range (cm)
HSL500/HSL510	X Small	40-50	28-32
	Small	50-60	30-36
	Medium	60-70	33-39
	Large	70-80	37-42
	X Large	80-90	40-46

SLING CODES

**Universal Deluxe Sling**

Product Code	Description	Material	Size
HSL100.PO.S	Universal Deluxe Sling with Loop	Polyester	Small
HSL100.PO.M	Universal Deluxe Sling with Loop	Polyester	Medium
HSL100.PO.L	Universal Deluxe Sling with Loop	Polyester	Large
HSL100.PO.XL	Universal Deluxe Sling with Loop	Polyester	Extra Large
HSL100.SS.S	Universal Deluxe Sling with Loop	Super Soft	Small
HSL100.SS.M	Universal Deluxe Sling with Loop	Super Soft	Medium
HSL100.SS.L	Universal Deluxe Sling with Loop	Super Soft	Large
HSL100.SS.XL	Universal Deluxe Sling with Loop	Super Soft	Extra Large
HSL110.PO.S	Universal Deluxe Sling with Clip	Polyester	Small
HSL110.PO.M	Universal Deluxe Sling with Clip	Polyester	Medium
HSL110.PO.L	Universal Deluxe Sling with Clip	Polyester	Large
HSL110.PO.XL	Universal Deluxe Sling with Clip	Polyester	Extra Large
HSL110.SS.S	Universal Deluxe Sling with Clip	Super Soft	Small
HSL110.SS.M	Universal Deluxe Sling with Clip	Super Soft	Medium
HSL110.SS.L	Universal Deluxe Sling with Clip	Super Soft	Large
HSL110.SS.XL	Universal Deluxe Sling with Clip	Super Soft	Extra Large
HSL300.PO.S	Universal Deluxe Sling with Head Support & Loop	Polyester	Small
HSL300.PO.M	Universal Deluxe Sling with Head Support & Loop	Polyester	Medium
HSL300.PO.L	Universal Deluxe Sling with Head Support & Loop	Polyester	Large
HSL300.PO.XL	Universal Deluxe Sling with Head Support & Loop	Polyester	Extra Large
HSL300.SS.S	Universal Deluxe Sling with Head Support & Loop	Super Soft	Small
HSL300.SS.M	Universal Deluxe Sling with Head Support & Loop	Super Soft	Medium
HSL300.SS.L	Universal Deluxe Sling with Head Support & Loop	Super Soft	Large
HSL300.SS.XL	Universal Deluxe Sling with Head Support & Loop	Super Soft	Extra Large

SLING CODES



In Chair Hammock Sling

Product Code	Description	Material	Size
HSL400.PO.L	In Chair Hammock Sling with Loop	Polyester	Small
HSL400.PO.L	In Chair Hammock Sling with Loop	Polyester	Medium
HSL400.PO.L	In Chair Hammock Sling with Loop	Polyester	Large
HSL400.PO.XL	In Chair Hammock Sling with Loop	Polyester	Extra Large
HSL400.SS.L	In Chair Hammock Sling with Loop	Super Soft	Small
HSL400.SS.L	In Chair Hammock Sling with Loop	Super Soft	Medium
HSL400.SS.L	In Chair Hammock Sling with Loop	Super Soft	Large
HSL400.SS.XL	In Chair Hammock Sling with Loop	Super Soft	Extra Large
HSL410.PO.S	In Chair Hammock Sling with Clip	Polyester	Small
HSL410.PO.M	In Chair Hammock Sling with Clip	Polyester	Medium
HSL410.PO.L	In Chair Hammock Sling with Clip	Polyester	Large
HSL410.PO.XL	In Chair Hammock Sling with Clip	Polyester	Extra Large
HSL410.SS.S	In Chair Hammock Sling with Clip	Super Soft	Small
HSL410.SS.M	In Chair Hammock Sling with Clip	Super Soft	Medium
HSL410.SS.L	In Chair Hammock Sling with Clip	Super Soft	Large
HSL410.SS.XL	In Chair Hammock Sling with Clip	Super Soft	Extra Large

SLING CODES



Toilet/Dress Sling

Product Code	Description	Material	Size
HSL200.PO.S	Toilet/Dress Sling with Loop	Polyester	Small
HSL200.PO.M	Toilet/Dress Sling with Loop	Polyester	Medium
HSL200.PO.L	Toilet/Dress Sling with Loop	Polyester	Large
HSL200.PO.XL	Toilet/Dress Sling with Loop	Polyester	Extra Large

Stand Aid Sling

Product Code	Description	Material	Size
HSL500.PO.S	Stand Aid Sling with Loop	Polyester	Small
HSL500.PO.M	Stand Aid Sling with Loop	Polyester	Medium
HSL500.PO.L	Stand Aid Sling with Loop	Polyester	Large
HSL500.PO.XL	Stand Aid Sling with Loop	Polyester	Extra Large
HSL510.PO.S	Stand Aid Sling with Clip	Polyester	Small
HSL510.PO.M	Stand Aid Sling with Clip	Polyester	Medium
HSL510.PO.L	Stand Aid Sling with Clip	Polyester	Large
HSL510.PO.XL	Stand Aid Sling with Clip	Polyester	Extra Large

MACKWORTH ESSENTIALS SATIN SHEET



The Mackworth Essentials Satin Sheet range is designed to reduce lifting associated strains and injuries to caregivers, whilst providing a solution for patients who require frequent repositioning.

The system consists of two components; a synthetic, breathable satin base sheet with a low friction central section and a soft top sheet with low friction backing. The Mackworth Essentials Satin Sheet range features handles on the Top Sheet to facilitate easier movement and use in conjunction with hoists.

KEY FEATURES & BENEFITS

Base Sheet

- Available in standard and bariatric mattress sizes
- Made of synthetic breathable satin and cotton

Top Sheet

- Maxi (full bed) and Midi (half bed) sizes
- Handles for easier positioning and use with hoists
- Made of synthetic breathable satin and cotton



Technical Information

Product Code	Description	Dimensions to fit mattress size
2002MAC9020	Mackworth Standard Base Sheet	2000(L)x900(W)mm
2002MAC1420	Mackworth Midi Top Sheet	2000(L)x1400(W)mm
2002MAC12020	Mackworth Satin Sheets Bariatric Base	2000(L)x1200(W)mm
2002MAC2020	Mackworth Satin Sheets Maxi Top	2000(L)x2000(W)mm

FREEWAY BED MANAGEMENT SYSTEM

FREEWAY BODY SUPPORT WEDGE

FREEWAY

The **Freeway Bed Management System** (bariatric) provides an easy to use, low friction solution for patients who require frequent repositioning. It consists of two separate units – a nylon, water resistant under sheet with wide non-slip sides, and a fully breathable, water resistant, quilted polyester softshell top sheet. The system has been designed to provide comfort for patients, and can help to reduce incidences of pressure ulcers. It also provides ease of use by carers, and forms an ideal bed solution as part of a single handed care package.

Also available to work in conjunction with the Freeway Bed Management System, is the **Freeway Support Wedge**. This provides a firm support for a patient in bed and helps them remain in a required position.

KEY FEATURES & BENEFITS

Bed Management System

- Provides comfort for patients
- Minimal effort from carers
- Ideal for use by a single carer
- Works easily with four section profiling beds
- Can be used on a pressure care mattress
- Can also be used with a hoist to assist turning

Body Support Wedge

- Provides a positive non-slip action when used against cotton bed linen
- Can be used to help reduce the incidence of pressure ulcers
- Made from 100% soft and flame retardant foam
- The polyurethane cover can be machine washed up to 80°C
- The nylon cover can be wiped clean with disinfectant



Technical Information

Product Code	Product Description	Dimensions
2002ACCBMS-COM	Freeway Bed Management System - Complete	-
2002ACCBMS-BSS	Freeway Bed Management System - Base Sheet Only	
2002ACCBMS-QS	Freeway Bed Management System - Top Sheet Only	2000(L)x1400(W)mm
2002ACCBMS-BCOM	Freeway Bed Management System - Complete - Plus Size	-
2002ACCBMS-BSB	Freeway Bed Management System - Base Sheet Only - Plus Size	2000(L)x1200(W)mm
2002ACCBMS-QB	Freeway Bed Management System - Top Sheet only	2000(L)x1900(W)mm
200FW-WEB53	Freeway Positioning Wedge	53mm(W)
200FW-WEB60B	Freeway Positioning Wedge - Plus Size	60mm(W)

PRISM MULTIGLOVE

PRISM LIMB LIFTER



The **Prism Multiglove** enables carers to facilitate a variety of moving and handling applications with ease due to its low friction. It is a simple but effective piece of equipment, and can be used to bring a person to a sitting position in a profiling bed, as the anti-shear material helps to smooth and release the tension on the skin.



Technical Information

Product Code	Description
2002ACCPM-MGLOVE	Prism Multiglove (Pair)

The **Prism Limb Lifter** has been designed to support limbs, allowing a number of procedures to be completed, i.e. wound care, without the carer having to support the limb. The Prism Limb Lifter can be easily attached to a hoist, enabling the carer to perform procedures in the safest possible working environment, eliminating the possibility of strain or injury from supporting the limb themselves.

Coloured inset loops are featured to aid correct positioning of the Prism Limb Lifter when attaching to a hoist. Each Prism Limb Lifter is supplied with two removable washable covers. Additional washable covers can be purchased here.

The product is Machine washable up to 72°C.



Technical Information

Product Code	Description	Dimensions
2002ACCPLL150X1	Prism Limb Lifter	150mm(L)
2002ACCPLL150X5	Prism Limb Lifter	150mm(L) (box of 5)

PRISM ONE WAY GLIDE PRISM HANDLING BELT



The **Prism One Way Glide Sheets** help a patient maintain a good sitting position and allows repositioning on a chair or bed.

The one-way glide will only slide in the direction of the arrows which are clearly labelled on the sheets. It can be used on all surfaces and helps prevent a patient from sliding forward in a wheelchair/day chair/bed and aids the epositioning of a patient into a good sitting position.




Technical Information

Product Code	Description
2002ACCPM-OWGS14	Prism One Way Glide - Single 1400 x 430
2002ACCPM-OWGS70	Prism One Way Glide - Single 700 x 430
2002ACCPM-OWGT60	Prism One Way Glide - Tube 600 x 430
2002ACCPM-OWGT73	Prism One Way Glide - Tube 730 x 600

The **Prism Comfi Belt** helps carers perform assisted manoeuvres and transfers safely whilst minimizing the risk of injury to themselves or the user. It has been designed to provide high levels of comfort to the user whilst providing the best possible grip and handling positions for the carer. Available in three separate sizes, each of which is also adjustable to ensure the perfect fit for waist sizes from 22inch (558mm) to 55inch (1395mm) .

The product is machine washable up to 72°C.



Technical Information

Product Code	Description
2002ACCPM-CBLARG	Prism Comfi Belt - Large
2002ACCPM-CBMED	Prism Comfi Belt - Medium
2002ACCPM-CBSMAL	Prism Comfi Belt - Small

TILT IN SPACE SHOWER CHAIR



This shower chair provides comfort and dignity for patients during washing and toileting routines, while maintaining ease of use for the carer.

It provides a versatile bathing and toileting solution for a wide variety of patients.

KEY FEATURES & BENEFITS

- Tilt-in-space shower/commode chair
- Adjustable seat angle of 6° to 25° ensures suitability for a wide range of conditions and facilitates patient transfers
- Length adjustable, swing-away footrest for ease of access and for varying leg lengths
- Height adjustable headrest suitable for a range of patient heights
- Flip-up armrests manufactured in strong yet lightweight aluminium
- Durable and easy to transport; will not rust or corrode in wet environments
- Padded seat provides a comfortable and secure seating position for the patient
- 4 individually locking castors provides a stable seating position whilst aiding manoeuvrability
- Removable backrest can be laundered and replaced as necessary

1

year warranty



158kg/25st
max user weight



Dimensions

Overall Width	650mm	Seat Width	450mm	Back Height From Seat	1000mm
Overall Height	1100-1400mm	Seat Height	550mm	Armrest Height From Seat	250mm
Overall Length	1000mm	Seat Depth	420mm	Armrest Width x Armrest Length	50 x 240mm

Technical Information

Product Code	Description
HMC119	Aluminium Tilt Transit Shower / Commode Chair

MACKWORTH ESSENTIALS SHOWER COMMODE CHAIR



The Mackworth Essentials shower chair is an attendant-propelled chair with height adjustable aluminium frame and a removable seat section to convert it from a moving, handling and showering solution to a commode chair. There is also an optional rear wheel and brake kit to simply and quickly transform it to a self-propelled chair. The chair benefits from the same rigorous, ergonomic design procedures and quality checks as other trusted Mackworth products. It is constructed from high quality materials that will provide years of reliable service.

KEY FEATURES & BENEFITS

- Cost effective solution without compromising on quality
- Height adjustable to fit all standard toilets
- Robust aluminium powder coated frame
- Removable back
- Self-propelled rear wheel kit available
- Swinging arms and foot-plates
- Removable commode access panel
- Round commode pan with handle fitted as standard
- Braked castors

1 year warranty

 130kg/20st
max user weight



Dimensions/Technical Information

Seat Dimensions	435(W) x 40(H) x 420(D)mm
Overall Width (attendant/self propelled)	562/620 mm
Chair Height (low, medium and high level settings)	965mm, 990mm, 1015mm
Back Height from Seat	440mm
Arm Height from Seat	235mm



Price Information

Product Code	Description
1601ESC0130	Mackworth Essentials Shower Chair

FREEWAY ADJUSTABLE HEIGHT SHOWER TROLLEY

FREEWAY

The Freeway Height Adjustable Shower Trolley benefits from optional cot sides, and its height can be adjusted by hydraulic foot pump. Available in different sizes of 1400mm, 1500mm, 1700mm and 1900mm, the Freeway Height Adjustable Shower Trolley has a working height of between 600mm and 925mm. The protective coating of all our shower trolleys and stretchers includes Freecote? ? our unique antimicrobial and antibacterial coating, which is resistant against the hospital super bug MRSA, e-coli, salmonella and listeria.

KEY FEATURES & BENEFITS

- Available in 3 sizes – 1500 mm, 1700 mm and 1900 mm
- Has a working height of 600 mm to 925 mm.
- Practical ergonomic design
- Durable, manoeuvrable and extraordinary value for money.
- Top quality stainless steel bearing castors
- Special custom build available
- Includes FREECOTE® – our unique antimicrobial and antibacterial coating

1 year warranty



200kg/30st
max user weight



With Adjustable
Backrest



Padded Changing
Trolley



With Liner

Technical Information

Product Code	Product Description
1801AHT1900	Adjustable Height Trolley - 1900mm
1801AHT1700	Adjustable Height Trolley - 1700mm
1801AHT1500	Adjustable Height Trolley - 1500mm
1801AHT1900B	Adjustable Height Trolley with Adjustable Backrest - 1900mm
1801AHT1700B	Adjustable Height Trolley with Adjustable Backrest - 1700mm
1801AHT1500B	Adjustable Height Trolley with Adjustable Backrest - 1500mm
1801AHTPC1900	Adjustable Height Padded Changing Trolley - 1900mm
1801AHTPC1700	Adjustable Height Padded Changing Trolley - 1700mm
1801AHTPC1500	Adjustable Height Padded Changing Trolley - 1500mm
1801AHT1900L	Adjustable Height Trolley with Liner - 1900mm
1801AHT1700L	Adjustable Height Trolley with Liner - 1700mm
1801AHT1500L	Adjustable Height Trolley with Liner - 1500mm

OPTIONAL EXTRAS:

- Low Cotside in 3 different heights (1300mm; 1802ACCWCU13, 1500mm; 1802ACCWCU15 and 1700mm; 1802ACCWCU17)
- Padding for Cotside; 1802ACCWCP

FREEWAY ASSISTED OR SELF-PROPELLED SHOWER, TOILET AND COMMODE CHAIRS

FREEWAY

Manufactured in the UK from only the highest quality materials, the comprehensive range of Freeway assisted or self-propelled shower, toilet and commode chairs are pioneering and firmly established as having a reputation second to none.

High quality, robust and durable, they are designed to provide the most effective solution to a client's toileting and showering needs.

Freeway Shower Chairs are of modular design, which means they can be configured to a choice of interchangeable seats, backrests, armrests and footrests.

A multitude of accessories are available to meet individual requirements and ensure the most bespoke solution.

All adult models within the range are specifically designed to fit over a standard toilet or commode, with three widths available:

A – 540 mm (21")

B – 490 mm (19")



C – 440 mm (17")

Paediatric, plus size and models to fit over automatic toilets are also included in the range.

The interchangeable makeup of Freeway Shower Chairs makes the range the most adaptable solution on the market, as accessories and parts can be added to suit the changing needs of a client.

Chairs can be manufactured with frames to be lowered and raised; rear wheels set back/forward and the frame depth extended and shortened.

Prism Medical UK has a team of experienced Assessors who can carry out a comprehensive shower chair assessment to determine the needs of the individual. To ensure the most effective and supportive solution to meet your client's needs, we highly recommend a shower chair assessment is made prior to purchase.



FREEWAY T80 RECLINING SHOWER CHAIR

FREEWAY

The Freeway T80 Shower Chair is highly versatile and the reclined position has significant advantages for increasing comfort and relieving pressure for a wide range of conditions, including:

- Multiple Sclerosis
- Spinal injuries
- Motor Neurone Disease
- Stroke & Other neurological conditions
- Cerebral Palsy

It has two variations; Freeway T80 COM and Freeway T80 AUTO.

Freeway T80 COM: A variation of the T80 reclining shower chair, designed to work seamlessly with the Clos-o-Mat Samoa toilet.

T80 AUTO: A version of the standard T80 which has been modified for use with automatic toilets including the Geberit 8000 range and Clos-o-Mat Palma Vita.

KEY FEATURES & BENEFITS

- This chair reclines through an infinite number of positions up to 45 degrees and is easy to operate.
- The T80 is available in self-propelled versions with a choice of reclines and comes with a wide range of positioning accessories.
- Featuring two high quality stainless steel gas struts, the tilt action is easy and light to operate.
- Available as assistant-propelled or self-propelled
- One year warranty

1 year warranty



Freeway T80 COM



Freeway T80 AUTO

Technical Information

Product Code	Description
1601T80C	T80 Tilt-in-Space Shower Chair – 440mm wide
1601T80B	T80 Tilt-in-Space Shower Chair – 490mm wide
1601T80A	T80 Tilt-in-Space Shower Chair – 540mm wide

FREEWAY T70 SHOWER CHAIR

FREEWAY

The Freeway T70 Shower Chair is available as a shower, toilet or commode chair. One of the major benefits of the Freeway T70 is that it is fitted with large-braked rear wheels and push rims making this chair extremely easy to self-propel.

As with all Freeway Shower Chairs, only the finest materials have been used during manufacturing to ensure the shower chair is durable and long lasting.

The Freeway T70 offers the ultimate in hygiene, as the durable seat, backrest and armrests can be simply wiped down after use. The protective coating of the shower chair includes FREECOTE® – our unique antimicrobial and antibacterial coating which is resistant to the hospital superbug MRSA, e-coli, salmonella and listeria.

KEY FEATURES & BENEFITS

- Available as a shower, toilet or commode chair
- Versions available to fit over Clos-o-Mat toilets, automatic bidets and for those who require easy side access
- Perfect for those who require a chair that is extremely easy to self propel
- Puncture proof rear braked wheels with removable push rims
- Large diameter castors to the rear allow easy manoeuvring over uneven surfaces
- Frame constructed from heavy duty, durable tubing in order to provide years of rust free performance
- Pliable backrest moulds to the user while providing excellent support
- Durable seat, backrest and armrests can be easily wiped clean to ensure excellent levels of hygiene
- Coating of the shower chair frame includes FREECOTE®, our unique antibacterial and antimicrobial coating

1

year warranty



200kg/30st
max user weight



Technical Information

Product Code	Description
1601T70C	Self propelled Shower/Toilet/Commode Chair – 440mm wide
1601T70B	Self propelled Shower/Toilet/Commode Chair – 490mm wide
1601T70A	Self propelled Shower/Toilet/Commode Chair – 540mm wide

FREEWAY T40 ASSISTANT PROPELLED SHOWER CHAIR

FREEWAY

The Freeway T40 is a highly versatile assistant propelled shower, toilet or commode chair that can be raised/lowered, extended forward, rear castors setback for amputee use.

KEY FEATURES & BENEFITS

- A highly versatile shower chair which can be made bespoke for a variety of applications
- Versions available to fit over Clos-o-Mat toilets and other automatic bidets
- Heavy duty total stop castors and footrests fitted as standard
- Frame constructed from heavy duty, durable tubing in order to provide years of rust free performance
- Durable seat, backrest and armrests can be easily wiped clean to ensure excellent levels of hygiene
- Coating of the shower chair frame includes FREECOTE®, our unique antibacterial and antimicrobial coating.

1 year warranty

200kg/30st
max user weight



Technical Information

Product Code	Description
1601T40C	T40 Assistant Propelled Shower Chair – 440mm wide
1601T40B	T40 Assistant Propelled Shower Chair – 490mm wide
1601T40A	T40 Assistant Propelled Shower Chair – 540mm wide

ELEXO SHOWER CHAIR



harvesthealthcare®
a prismhealthcare company

The Elexo is a height-adjustable shower-toilet chair which provides comfort and security for both carers and patients. Built for safety, comfort and cleanliness, the Elexo is an exceptional tool for carers looking to give residents a secure and calm showering or toileting experience.

KEY FEATURES & BENEFITS

- Simple electric height adjustment
- Hand control with battery indicator
- Sleep mode setting after 60 minutes
- Tilt adjustment to a maximum of 17 degrees
- Easy clean, hygienic seat with toilet opening
- Inlay for the toilet opening
- Foldable upholstered armrests with front lock to enable sideways transfer
- Ergonomically formed seat and back rests
- Fold-out footrests with calf protection
- Fold-away and removable leg supports
- Rail for toilet bucket or bedpan
- Durable castors with swivel & lock features
- Extremely stable due to the fitting of two lifting motors
- Simple to recharge, with charger included

3 year warranty



180kg/28st
max user weight



**Ergonomically
formed seat for
toilet use**

Technical Information

Product Code	Seat Width	Seat Depth	Adjustable Height
HLA013	480mm	460mm	50m-105mm

AUTOLIFT BATHING AID



harvesthealthcare®
a prismhealthcare company

The Autolift is a strong, attractively designed bathing aid which utilises very simple principles to solve the bathing problems experienced by the elderly and disabled.

High quality materials are used in the construction of Autolifts, to help withstand humid bathroom atmospheres and give years of reliable service.

KEY FEATURES & BENEFITS

- Consists of three basic parts: a mounting plate that is mounted on the floor next to the bath, a mechanically operated lift column with arm and a comfortable, shaped chair
- Available in both electric and manual version
- Can be fixed either to the side or end of the bath
- Various arms make it suitable for almost any bathroom
- Can be used with either a fixed or a detachable seat option with additional sub-chassis
- By using solid and durable materials, the Autolift is suitable for years of service in wet environments

3 year warranty



127kg/20st
max user weight



OPTIONAL EXTRAS:

- Headrest
- Biocompatible safety belt
- Extra detachable battery pack and charger for the electric version

Technical Information

Product Code	Short Arm (Left to Right)	Long Arm (Left to Right)
HLA025	430mm	480mm

EVE

SHOWER AND BATH LIFTER



The sleek Eve Shower and Bath Lifter allows users to be able to safely and comfortably shower or bathe. Highly manoeuvrable, and fitting in most modern care bathtubs, the Eve is suitable for residents who need extra support getting in and out of the bath and is also fully functional as a shower chair.

KEY FEATURES & BENEFITS

- Made entirely of aluminium and stainless steel
- 100% functional as a shower chair
- Simple operation by only one carer
- Unrestricted transfer from the Ascent Pro 200 to the seat lifter
- Ergonomically designed safety steering handles ensure accurate manoeuvring
- Easy-to-maneuver double castors
- Handheld control unit with visual indication of battery level
- Integrated service module reports technical problems with the lift as well as necessary maintenance.

3

year warranty

185kg/29st
max user weight



OPTIONAL EXTRAS:

- Digital scale
- Replacement battery
- External wall charging station

Technical Information

Product Code	Length	Width	Height	Lifting Range	Free Space Under Carrier Frame
HLA008	950mm	680mm	1240mm	495mm-1065mm	30mm



harvesthealthcare®
a prismhealthcare company

Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

T +44 (0)1709 377172
F +44 (0)1709 377173
E sales@harvesthealthcare.co.uk

www.harvesthealthcare.co.uk

Plus Size Brochure - May 2021

© Copyright Harvest Healthcare

Harvest Healthcare Limited. Company No: 07210261

Our Full Terms & Conditions are available by request or can be found on our Website.

Harvest Healthcare reserves the right to alter or amend this document without notice.

