

Woburn Profiling Bed



Product Code:	HLB791.03
Product Name:	Woburn Profiling Bed
Description:	<p>Robust 115mm wooden side rails collapse below the mattress surface level when lowered, and extend higher to allow a safe mattress depth of 205mm, giving scope for improved pressure relief and comfort.</p> <p>The quality features and reliability continue to provide outstanding comfort and ease of care for occupant and carer.</p> <p>Adjustment of back rest, knee break and platform height all at the touch of a button.</p>
Safe Working Load:	220kg / 34 Stone
Max. Patient Weight:	185kg / 29 Stone
Overall Dimensions:	2135 x 1025mm
Mattress Base:	2000 x 900mm
Max. Mattress Depth:	205mm
Adjustable Height:	365 - 795mm
Angle of Back Rest:	Up to 70 degrees
Angle of Knee Break:	Up to 30 degrees
Bed Weight:	86kg + 10kg with rails
Key Information:	<p>4 Section Profiling with Optional Auto-Regression, Trendelenburg function, in-built extension and lifting pole</p> <p>Wraparound Wooden Head and Foot Boards</p> <p>Electronically Controlled Back Rest and Knee Break</p> <p>Integral Wooden Side Rails</p> <p>Manually Adjusted Leg Rest</p> <p>Low Voltage – 24v; 9v Battery Backup</p> <p>Lockable Handset</p> <p>Comprehensive manual; CE marked</p>

Supplier Contact Details:

t: +44(0)1709 377172

e: enquiries@harvesthealthcare.co.uk

w: www.harvesthealthcare.co.uk

IMPORTANT SAFETY NOTICE

HARVEST HEALTHCARE LTD SUPPLIES EQUIPMENT AND APPLIANCES TO COVER A RANGE OF MEDICAL CONDITIONS. SPECIFIC CHOICE OF EQUIPMENT IS THE SOLE RESPONSIBILITY OF THE CLINICAL PRACTITIONER OR CARER OVERSEE-ING INDIVIDUAL PATIENT CARE.

RISK ASSESSMENTS SPECIFIC TO USER MUST BE CARRIED OUT PRIOR TO THE USE OF ALL PRODUCTS. PRODUCTS AND DATA CONTINUALLY UNDERGO DEVELOPMENT AND MAY THEREFORE DEVIATE FROM THE DETAILS GIVEN.

ISSUE 2: DEC 2017