

WOBURN LOW PROFILING BED



Product Code:	HLB797.03
Description:	The dependable 4 section profiling bed now features high quality control system, a reduced minimum platform height and a higher load capacity. By lowering the platform to only 220mm from floor level, the risk of injury for those with a tendency to roll out of bed is significantly reduced. Comfort setting of height, back rest and knee break angles is easy when using the lockable handset.
Safe Working Load:	220kg (34 Stone)
Max. Patient Weight:	185kg (29 Stone)
Mattress Base Size:	2000 x 900mm
Max. Mattress Depth:	205mm
Adjustable Height:	220 - 650mm
Angle of Back Rest:	Up to 70 degrees
Angle of Knee Break:	Up to 30 degrees
Bed Weight:	89.5kg (without side rails)
Key Information:	<ul style="list-style-type: none"> 4 Section Profiling Bed Electronically Powered Back Rest & Thigh Rest Manual Leg Rest Adjustment Wraparound Wooden Head and Foot Boards Side Rails - Optional Trendelenburg Function and Auto-Regression - Optional In-Built Extension - Optional Lockable Handset Low Voltage – 29v; 9v Battery Backup Comprehensive Manual CE marked

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

Important Safety Notice

< Ufj Ygh< YUth VWfY @fX gi dd`Yg Yei Jda YbhUbX`Udd`JubWg hc`Wvj Yf`U`fub[Y`cZa YX]WU`VebX]hcbg" GdYVW]MWKc]Mf`cZyei Jda Ybhjg`h.Y gc`Y`fygdcbg]M]]mcZh.Y Wb]WU`dfUWf]c:bYf cf`WUfYf`
cj`YfgY`Y]b[]bX]]xi`U`dUf]bhVWfY" Rig_`UggYgga`Ybtg`gdYVW]Wc`i`gYfa`i`ghVY`WUff]YX`ci`hd]c`f`hc`h.Y`i`gY`cZU`dfcXi`WfY"
PfcXi`WfY`UbX`XUfU`Veb]bi`U`mi`bXYf[c`XYj`Y`cda`YbhUbX`a`Umh.YfY`Z`fY`XYj`]UfY`Zca`h.Y`XYfU]g[]`Yb"

Version 2: Dec. 2017