

WOBURN COMMUNITY PROFILING BED



Product Code:

HLBC791.02

Description:

- Robustly designed for the rigours of community use this 4 Section Profiling Bed builds on the quality, technologies and features of our renowned Woburn range - at a lower cost.
- Top quality components and high capacity.
- Full length rails.
- Ideal for many environments, the bed is easily cleaned, serviced and handled.
- Low voltage to power a range of accessories including mattresses.
- Extendable base

Safe Working Load:

220kg / 34 Stone

Max. Patient Weight:

185kg / 29 Stone

Overall Dimension:

2115 x 1015mm

Mattress Base:

2000 x 900mm

Max. Mattress Depth:

205mm (when side rails are fitted)

Adjustable Height:

385mm-810mm

Angle of Back Rest:

Up to 70 degrees

Angle of Knee Break:

Up to 30 degrees

Bed Weight:

81kg + 11kg with rails

Side Rails

Description:

Integral Wooden Side Rails with plastic end caps

Dimensions:

1973 x 115 x 28mm (with end caps)

Key Information:

Trendelenburg Function as standard

Auto- Regression function and Inbuilt Extension are available as options

Manual Leg Rest Adjustment

Low Voltage – 29v; with 9v Battery Backup

Comprehensive User manual

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60