

WOBURN COMMUNITY LOW PROFILING BED



Product Code:	HLBC791L HLBC791LAR (with Auto-Regression)
Description:	<ul style="list-style-type: none"> • Built to the same specification as our popular Woburn Community Bed • Offers the benefit of a lower entry and exit height • Ideal for many environments, fully compliant with legislation. • Comfort adjustments of height, back rest or knee break done simply using the lockable handset.
Safe Working Load:	220kg / 34 Stone
Max. Patient Weight:	185kg / 29 Stone
Overall Dimension:	2115 x 1015mm
Mattress Base:	2000 x 900mm
Max. Mattress Depth:	205mm
Adjustable Height:	220mm-640mm
Angle of Back Rest:	Up to 70 degrees
Angle of Knee Break:	Up to 30 degrees
Bed Weight:	81kg + 11kg with rails
Key Information:	4 Section Profiling Bed Electronically Powered Back Rest & Thigh Rest Manual Leg Rest Adjustment Wooden Head and Foot Boards Side Rails - Optional Auto-Regression and Inbuilt Length Extension - Optional Trendelenburg Function as standard Lockable Handset Low Voltage – 29v; 9v Battery Backup Comprehensive Manual Warranty: 5yrs

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk
 Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

Important Safety Notice

<Ufj Ygh< YUth.VWfY @fX gi dd]Yg Yei]da YbhUbX'Udd.]ubWg hc'Wj Yf'U'fub[Y'cZa YX]WU'W'ebX]hc'bg" GdYVW]MwK.c]M' cZYeI]da Ybhjg'h.Y gc'Y fYgdcbg]M]]mcZk.Y W]b]WU' d'fUW]j]c:byf cf'WUfYf
cj YfgYf]b[]bX]]xi U' dUf]bhVWfY" Rig_ UggYgga Ybtg'gdYVW]M'c: i gYfa i ghVY WUff]YX'ci hdf]cf hc'h.Y i gY'cZU' d'fcXi V]g"

Version 2: Nov 2017