

# SPECIALIST ANTI-ENTRAPMENT BUMPER FOR SMALL AND JUNIOR ADULTS



<b>Product Code:</b>	HJB002
<b>Description:</b>	<p>Specialist Anti Entrapment Bumper</p> <ul style="list-style-type: none"> <li>• The Specialist Anti Entrapment Bumper has been designed to fit onto a Woburn community bed.</li> <li>• The product fills the need to provide a safe environment for a junior or a vulnerable small adult who has an adult-sized bed but needs to meet the new safety legislation BS EN50637:2017 for junior beds.</li> <li>• This product conforms to BSEN50637:2017 and EN60601-2-52 ideal for long term care, such as children's nursing homes and care in a domestic environment. To meet this standard fully you must also purchase a separate locking box (product code: BOX001) to restrict unwanted access to the handset.</li> <li>• The Specialist Anti Entrapment bumper prevents access to finger trap zones and allows the occupant to see out, which in turns helps the carer as they can see through the netting without having to disturb the user.</li> </ul>
<b>Note:</b>	The Specialist Anti Entrapment system must not be used by a person smaller than 125cm.
<b>Cleaning and Care:</b>	<p>The cover is removable and machine washable, for better infection control. Pre wash cold 10 minutes Main wash 80°C 10 minutes Followed by cold rinses and extraction. Do not use abrasive cleaners, phenol disinfectants, solvents, or alcohol-based cleansers, e.g. Dettol, Phenicol, Hibiscrub, Clearsol, Stericol, Hycoline, as these will damage the cover materials. Tumble dry at low heat. Do not iron. Ensure that the cover is thoroughly dried before inserting the foam and placing the bumper into storage.</p>
<b>Latex free: User</b>	Yes
<b>Instruction:</b>	Yes
<b>Flammability:</b>	BS 7175: 1989 (2013) ignition sources 0,1

## Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk  
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

## Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 2: March 2022