

ULTIMATE 1200 LIFTING POLE



Product Code: HLDU582.12

- Description:**
- Optional lifting pole is compatible with the Woburn Ultimate 1200 Bed.
 - slots securely into sockets provided in the bed frame.
 - Has a durable powder coated finish.
 - Assists the occupant to sit up in bed by holding onto the triangular handle.

Safe Working Load: 80 kg

Weight: 11.9 kg

Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Lifting Pole is guaranteed for one year from date of purchase.

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 2: Nov 2017