Sling Assessment

Selecting the appropriate sling size for your resident

Name


Company:


Address:

Telephone:

Email:

HSL100/110 Universal Deluxe - Loops/Clips







All measurements are in cm.

Using the Sling Sizing Chart

Step 1:

Measure the width of your resident between the outer midpoints of each upper arm as per the diagram (W).

Step 2:

Input your measurement into the corresponding box under 'User Width', next to to the width range it fits within.

Step 3:

Measure the back length of your resident as per the diagram (B).

Step 4:

Input your measurement into the corresponding box under 'User Back Length', next to the width range it fits within. Please note some ranges overlap - in these cases select the most appropriate size.

Step 5:

If these two measurements correlate your resident fits a standard sling. If not, consider a bespoke sling or select the most appropriate standard sling after a risk assessment.

All sling selections should be accompanied by a risk assessment by a trained person. People of atypical anatomy may not be suited to a standard sling.

Material

Polyester

Super Soft

Comments / Additional info:

Date:



t +44 (0)1709 377172 f +44 (0)1709 377173 e sales@harvesthealthcare.co.uk www.harvesthealthcare.co.uk

Sling Sizing Chart

48-56 53-63 63-73 68-78 68-78
63-73
68-78
68-78
Clip Sling

Signature: