

# HEAD AND FOOTBOARD BUMPER



## Product Code:

**HBCWOB:** Head and Footboard Bumper for Woburn Bed

**HBCCOM1.12:** Head and Footboard Bumper for Woburn Community 1200 Bed

**HBCULT:** Head and Footboard Bumper for Woburn Ultimate Bed

## Description:

- The Head and Footboard bumper fits all Woburn and Community Beds.
- Providing extra protection for the user at the head and footend of the bed.

## Cleaning and Care:

The cover is removable and machine washable, for better infection control. Pre wash cold 10 minutes

Main wash 80°C 10 minutes

Followed by cold rinses and extraction.

Do not use abrasive cleaners, phenol disinfectants, solvents, or alcohol-based cleansers, e.g. Dettol, Phenicol, Hibiscrub, Clearsol, Stericol, Hycoline, as these will damage the cover materials.

Tumble dry at low heat.

Do not iron.

Ensure that the cover is thoroughly dried before inserting the foam and placing the bumper into storage.

## Latex free:

Yes

## Flammability:

BS 7175: 1989 (2013) ignition sources 0,1

## Supplier Contact Details:

t: +44(0)1709 377172

e: [enquiries@harvesthealthcare.co.uk](mailto:enquiries@harvesthealthcare.co.uk)

w: [www.harvesthealthcare.co.uk](http://www.harvesthealthcare.co.uk)

Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

## Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products.

Products and data continually undergo development and may therefore deviate from the details given.

Version 1: March 2022