

WOBURN ULTRA-LOW 800 PROFILING BED



Product Code:	HLB800.01AR
Description:	<p>A fully functioning four- section Ultra-Low profiling bed. Platform can be lowered to 90mm off the floor for patient safety and raised to 800mm to make it easier for the carer when providing care for the patient.</p> <p>Electronically setting of height, back rest or knee brake is simple, using the lockable handset.</p> <p>The Woburn Ultra-Low 800 has a Trendelenburg function as standard. Optional bed levels, lighting and accessories available.</p>
Max Patient Weight:	175kg (27 Stone)
Safe Working Load:	185kg (29 Stone)
Overall Dimension:	2399 x 900mm
Adjustable Platform Height:	90mm -800mm
Mattress Dimensions:	2000 x 900mm
Angle of Back Rest:	Up to 70 degrees
Angle of Knee Brake	Up to 30 degrees
Bed Weight:	125 kg
Key Information:	<ul style="list-style-type: none"> • Low voltage • Lockable handset • Emergency battery backup • Wooden head and foot boards • In-built auto-regression • Trendelenburg as standard • Electric backrest and kneebreak with mechanical lower leg support • Easy Cleaning - complies to MHRA guidelines • Compatible with environmental control system • Comes with a 5-year warranty

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

Important Safety Notice

<Ufj Ygh< YUth VUfY @fX gi dd`Yg Yei Jda YbhUbX Udd`JubWg hc`Wj Yf U flub[Y`cZa YX]WU`VebX]hcbg" GdYVW]Wk.c]Mf`cZYeI Jda Ybhig`h.Y gc`Y fYgdcbg]M]]mcZh.Y Wb]WU`dfUW]j:byf cf WUfYf
cj YfgY]b[]bX]]xi U`dU]YbhVUfY" Rig_ UggYgga Ybtg`gdYVW]Wk`i gYfa i ghVY WUff]YX`ci hdf]cf hc`h.Y i gY`cZU`dfcXi V]g`
PfcXi V]g UbX XUfU Veb]bi U`mi bXYf[c XYj Y cda YbhUbX`a Umh.YfY`Z:fY XYj]UfY Zca` h.Y XYU]g [] Yb`

Version 3: May 2022