

AirSet Cushion

**Product Code:****HHC001****Description:**

The Harvest Healthcare AirSet Cushion is a state-of-the-art pressure care solution that boasts an innovative hybrid design, combining the benefits of foam and air to provide optimal pressure relief.

Users can easily adjust the cushion to their desired level of pressure relief, thanks to its self-inflating technology, ensuring maximum comfort and support. The cushion is climate-controlled, helping maintain a comfortable seating temperature, reducing the risk of heat build-up and pressure ulcers. The AirSet Cushion's unique construction allows it to distribute pressure evenly across the seating surface, promoting healthy circulation and reducing the risk of injury.

| | |
|------------------------------|----------------------------|
| Dimensions: | 40cm(L) X 40cm(W) X 7cm(H) |
| Weight of Cushion: | 1.5kg/0.2st. |
| Max. User Weight: | 120kg/18.8st. |
| Pressure Relief Risk: | High-Risk |
| Number of Cells: | 36 |

Warranty

In terms of Harvest Healthcare Guarantee and Warranties the AirSet Cushion are guaranteed for three years from date of purchase.

Exclusions:

1. Removable cushion covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded
3. Components of pumps which are recommended for annual service, are guaranteed for one year.

Supplier Contact Details:

t: +44(0)1709 377172

e: sales@harvesthealthcare.co.uk

w: www.harvesthealthcare.co.uk

Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact.

(Mattress risk categories are to be used as a guide only)

Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 1: April 2023

