



# VANDAL PROOF FIRE-TECH 7 MATTRESS



**Product Codes:** HSMFIRET.VP

- Description:**
- The Fire-Tech 7 is a cost-effective high quality mattress designed for challenging environments where you need a mattress that is fully secure and fire resistant.
  - Has a Crib 7 fire rating - the most demanding fire test in the industry.
  - The mattress comes complete with a welded, anti-ligature, water resistant and fully-wipable cover.

**Mattress Dimensions:** 2000 x 900 x 150mm  
1900 x 900 x 150mm

**Mattress Weight:** 8.5kg

**Max. User Weight:** 30 Stone (190kg)

**Density Range:** 31-34 (kg/m<sup>3</sup>)

**Foam Hardness Range:** 100-130N

**Fire Retardancy:** BS 7177 Ignition Source 7

**Cleaning Guidance**

The mattress foam may be hand washed with mild detergent. Do not immerse in water. Allow to dry naturally. Cover should be wiped down with hot water and mild detergent or machine washed up to 95°C and tumble dried. For a comprehensive list of tested cleaning products see our website.

**Warranty:**

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Fire-Tech 7 Vandal Proof mattress is guaranteed for three years from date of purchase.

**Exclusions:**

1. Removable covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded

**Supplier Contact Details:**

t: +44(0)1709 377172 e: [enquiries@harvesthealthcare.co.uk](mailto:enquiries@harvesthealthcare.co.uk); w: [www.harvesthealthcare.co.uk](http://www.harvesthealthcare.co.uk)  
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

**Important Safety Notice**

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact.  
(Mattress risk categories are to be used as a guide only) Risk assessments specific to user must be carried out prior to the use of all products.  
Products and data continually undergo development and may therefore deviate from the details given.