

REFLECT 2 MATTRESS



Product Codes:

HFR147
HFR147/12
HFR147/1370

Description:

- Castellated viscoelastic topped and firm based profiling foam mattress.
- Castellated viscoelastic foam for improved immersion and pressure distribution.
- Firm base for support and evenly distribution of weight.
- Keyhole slots in the base aid profiling.
- Four-way stretch, water resistant, anti bacterial zipped cover with welded seams.
- No turn mattress.
- Available in 3 mattress widths.

Mattress Dimensions and Weight, per size:

2000 x 900 x 165mm	15kg
2000 x 1200 x 165m	22kg
2000 x 1370 x 165mm	25kg

Max. User Weight:

39 Stone (247kg)

Risk Level:

High

Cover:
(Available Separately)

Four-way stretch, welded, vapour permeable, PU coated cover. Silver-based antibacterial properties in the PU coating are washable. Cover does not contain latex.

Cleaning Guidance

The cushion foam may be hand washed with mild detergent. Do not immerse in water. Allow to dry naturally. Cover should be wiped down with hot water and mild detergent or machine washed up to 95°C and tumble dried. For a comprehensive list of tested cleaning products see our website.

Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Reflect 2 mattress is guaranteed for three years from date of purchase.

Exclusions:

1. Removable covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact.

(Mattress risk categories are to be used as a guide only) Risk assessments specific to user must be carried out prior to the use of all products.

Products and data continually undergo development and may therefore deviate from the details given.