


COMMUNITY MATTRESS



Product Codes:

HSMCOM
HSMCOM.12

Description:

- The Community Mattress is a cost effective comfort static mattress for use on profiling beds.
- The high resilient, fire retardant foam is covered with a zipped four-way stretch, water resistant PU cover and is shaped to allow profiling.
- The community mattress provides excellent pressure distribution and can be turned to extend service life.

Mattress Dimensions:

2000 x 900 x 125mm
2000 x 1200 x 125mm

Mattress Weight:

8kg
10kg

Max. User Weight: Risk

25 Stone (159kg)

Level:

Medium

Foam Density Range:

31-34 (kg/m³)

Foam Hardness Range:

100-130N

Fire Retardancy:

BS 7177 Ignition Source 5

Cover:

(Also available separately)

Four-way stretch PU coated polyester cover, with zip fastening. The coating is water resistant and antibacterial as well as being vapour permeable, washable and flame retardant (crib 5).

Cleaning Guidance

The mattress foam may be hand washed with mild detergent. Do not immerse in water. Allow to dry naturally. Cover should be wiped down with hot water and mild detergent or machine washed up to 95°C and tumble dried. For a comprehensive list of tested cleaning products see our website.

Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Community mattress is guaranteed for three years from date of purchase.

Exclusions:

1. Removable covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact.

(Mattress risk categories are to be used as a guide only) Risk assessments specific to user must be carried out prior to the use of all products.

Products and data continually undergo development and may therefore deviate from the details given.