


COMMUNITY PLUS MATTRESS



Product Codes:

HSMCOMP
HSMCOMP.12
HSMCOMP.DV

Description:

- The Community Plus Mattress is a higher quality mattress but remains very cost-effective
- The mattress foam is thicker than the Community mattress and is highly resilient, providing better pressure distribution and extra comfort
- Comes complete with a four-way stretch, antibacterial, water resistant and washable PU coated cover.
- Additional sizes available on request

Mattress Dimensions:

2000 x 900 x 150mm-HSMCOMP
2000 x 1200 x 150mm-HSMCOMP.12
1900 x 900 x 150mm-HSMCOMP.DV

Mattress Weight:

8.5kg-HSMCOMP & HSMCOMP.DV
11kg-HSMCOMP.12

Max. User Weight:

30 Stone (190kg)

Risk Level:

High

Waterlow Score:

15-19

Foam Density Range:

31-34 (kg/m³)

Foam Hardness Range:

100-130N

Fire Retardancy:

BS 7177 Ignition Source 5

Cover:

(Also available separately)

Four-way stretch PU coated polyester cover, with zip fastening. The coating is water resistant and antibacterial as well as being vapour permeable, washable and flame retardant (crib 5).

To be used in conjunction with 2 hourly turn

Cleaning Guidance

The mattress foam may be hand washed with mild detergent. Do not immerse in water. Allow to dry naturally. Cover should be wiped down with hot water and mild detergent or machine washed up to 95°C and tumble dried. For a comprehensive list of tested cleaning products see our website.

Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Community Plus mattress is guaranteed for three years from date of purchase.

Exclusions:

1. Removable covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact.

(Mattress risk categories are to be used as a guide only) Risk assessments specific to user must be carried out prior to the use of all products.

Products and data continually undergo development and may therefore deviate from the details given.