



harvesthealthcare[®]
a prismhealthcare company

Cadence

Bed Range



**BROWSE
THE RANGE**

www.harvesthealthcare.com





INTRODUCING **Cadence**

Get ready to revolutionise your care environment with Harvest Healthcare's newest addition - the Cadence Bed Range!

With cutting-edge technology and comfort at the forefront, our Cadence Beds are the ultimate solution for resident care.

Experience unparalleled comfort and style with our advanced features and state-of-the-art technology. But that's not all - the Cadence Bed Range is also fully customisable to your exact specifications. Choose from a wide range of bold colours and calming fabrics to make your bed the perfect addition to any care environment.

Design the perfect bed and elevate your resident's well-being and safety with the Harvest Healthcare Cadence Bed Range.

Don't wait - create your bespoke solution today!

TABLE OF CONTENTS

Bed Configuration	1
Bed Width	2
Bed Frame	3
Headboard & Footboard	4
Colours & Fabrics	6
Side Rails	8
Side Skirts	9
Pre-Configured Range	10
Cadence Comfort	11
Cadence Select	12
Cadence Premium	13
Bed Specifications	14
Why Choose Harvest?	16
Contact Us	See Back

BED CONFIGURATION

Design Your Perfect Bed



Cadence Configuration

BED WIDTH

Sleep Comfortably: Choose Your
Desired Bed Width.



900mm
SYMA301



1050mm
SYMA302

Cadence Configuration

BED FRAME

Sleep Soundly: Choose your ideal bed frame.



Standard Frame
SYMA301



With Central Brake



Cadence Configuration

HEADBOARD & FOOTBOARD


Customise Your Head & Footboard:

The Cadence Bed Range offers a variety of headboard and footboard options to choose from, so you can create a bed that's both stylish and functional. Whether you're looking for a traditional headboard with a clean finish or a more modern footboard with a sleek design, we have something to fit your taste.


Our headboards and footboards are made from high-quality materials that are built to last. They're also easy to care for, so you can keep your bed looking its best for years to come.

See next page for headboard options


900mm/1050mm



Standard Headboard C
SYMB309/SYMB315



Standard Headboard V
SYMB307/SYMB311



Curved Headboard
SYMB301/SYMB313



Low Footboard V
SYMB308/SYMB312



Curved Low Footboard
SYMB305/SYMB314



Cadence Configuration

COLOURS & FABRICS

Design your perfect look. Choose from a wide range of colours, fabrics, and upholstery options to create a unique and personalised style for your space. Mix and match to find the perfect combination that suits your taste and complements your decor.

We also offer the option to colour-match your new bed to your existing furniture.

Please get in touch for upholstery options



Walnut



Truffle Sonoma Oak



Oyster Urban Oak



Mountain Oak



Grey Nebraska Oak



Grey Bardolino Oak



Lindberg Oak



Lissa Oak

Cadence Configuration

SIDE RAILS

Sleep Securely: Choose your ideal side rails.



2-Part Split Siderails
SYMD304



3-Part Split Siderails
SYMD303



Harvest Siderails
HLDWC110

Cadence Configuration

SIDE SKIRTS

Choose from a variety of styles, colours, and materials to create a look that perfectly matches your decor. Whether you prefer a classic or modern aesthetic, we have the perfect side skirt to complement your bed.

Side skirts are available as an optional extra.



Standard Side Skirt
SYMC301

PRE-CONFIGURED RANGE

Orchestrated to Perfection






Cadence Comfort


CADENCE.COMFORT.C

The Cadence Comfort Bed is a stylish, versatile and cost-effective option that features a standard head and footboard with a curved design, complete with side skirts. The bed can be used with either continuous side rails or split side rails, which are available as optional extras, providing the flexibility to meet a variety of needs and preferences.


With its elegant design and adaptable features, the Cadence Comfort Bed is a comfortable and stylish choice for any bedroom.



Standard Headboard C or Standard Headboard V




(Footboard will be the same as selected Headboard)



Standard Frame



Harvest Siderails (Optional)



Standard Sideskirts (Optional)

FOR FULL SPECIFICATIONS PLEASE SEE PAGE 14




Cadence Select


CADENCE.SELECT.LFB

The Cadence Select bed is a sleek and functional mid-range option that boasts a standard headboard and a low footboard, both adorned with side skirts for added style. The bed is designed to be used with split side rails, which are available as optional extras, providing additional support and customisation options to suit individual needs.


With its sophisticated design and versatile features, the Cadence Select bed is an ideal choice for those looking for a comfortable and practical bed frame that doesn't sacrifice style.



Standard Headboard V




Low Footboard V



Standard Frame



3-Part Split Siderails (Optional)



Standard Sideskirts (Optional)

FOR FULL SPECIFICATIONS PLEASE SEE PAGE 14








Cadence Premium

CADENCE.PREMIUM

The Cadence Premium bed is a luxurious and elegant bed frame that boasts a curved headboard and low curved footboard, both accented with stylish side skirts. The bed is designed to be used with split side rails, which can be added as optional extras, providing enhanced support and customisation options to meet individual needs.

With its refined design and customisable features, the Cadence Premium bed is a perfect choice for those seeking a high-end and personalised sleeping experience.

-  → Curved Headboard
-  → Curved Low Footboard
-  → Standard Frame
-  → 2-Part Split Siderails
(Alternative siderails are available)
-  → 2-Part Split Sideskirts
(Alternative sideskirts are available)

FOR FULL SPECIFICATIONS PLEASE SEE PAGE 14











Cadence Specifications

The “all-in-one care bed” combines all the requirements of a modern care bed in terms of comfort, ergonomics and aesthetics.

This section of the brochure highlights the key specifications of the Cadence Bed range, including the bed's dimensions, weight capacity, and mobility features.

It also provides information about the bed's advanced positioning capabilities, including Trendelenburg and auto-regression, as well as the ability to adjust the height and angle of the bed to suit the needs of the user.

Key Features

	Warranty	5 Years
	Trendelenburg	Yes
	Auto-Regression	Yes
	Max User Weight	200kg
	Safe Working Load	235kg
	Dimensions (mm)	2155 x 1055
	Height Range (mm)	240 - 800
	Bed Weight	126kg

FOR MORE INFORMATION VISIT OUR WEBSITE



Cadence Specs Cont.

A comprehensive overview of the advanced features and functions that set our Cadence beds apart from the competition.



Backrest
0 - 70°



Footrest
0 - 45°



Trendelenburg
Included



Height Range
240 - 800mm

Additional Features

- ✓ Extendable
- ✓ Standard & Extended Side Rails
- ✓ 70 Degree Back Rest Angle
- ✓ 45 Degree Thigh Rest Angle
- ✓ Lockable Handset
- ✓ Low Voltage Power Supply
- ✓ Emergency Battery Back-Up
- ✓ 5 Position Foot Rest
- ✓ IPX4 Water Resistance
- ✓ Reverse Trendelenburg

FOR MORE INFORMATION VISIT OUR WEBSITE

WHY CHOOSE HARVEST?



Upgrade or change the look and feel without buying a new bed frame.



In-house upholstery and design team offers unrivalled design flexibility.



Our products are of exceptional quality and come with a 5-year warranty.




Nationwide engineering team provides lifetime servicing and asset register provision.



MO training available for bed care and troubleshooting.

DELIVERING COMFORT, CARE AND CONVENIENCE WITH EVERY ORDER.



Commitment to honest and friendly service.

Responsive and efficient in ensuring your needs are met.

Delivery team goes the extra mile, providing installation and setup services on a date of your choosing.

Maintenance backup service, with readily available parts for optimal performance over time.



harvesthealthcare®
a prismhealthcare company

GET IN TOUCH



+44 (0)1709 377172



www.harvesthealthcare.co.uk



sales@harvesthealthcare.co.uk



Sheaf House, Bradmarsh
Business Park, Bradmarsh
Way, Rotherham S60 1BW



[harvest-healthcare-ltd](https://www.linkedin.com/company/harvest-healthcare-ltd)



[@HarvestHealthUK](https://twitter.com/HarvestHealthUK)



[Harvest Healthcare](https://www.youtube.com/channel/UC...)



Visit the
Cadence
Webpage

Harvest Healthcare Cadence Bed Range - April 2023
© Copyright Harvest Healthcare
Harvest Healthcare Limited. Company No: 07210261

Our Full Terms & Conditions are available by request or can be found on our Website.
Harvest Healthcare reserves the right to alter or amend this document without notice.