

Cadence Comfort Profiling Bed



Product Code:

CADENCE.COMFORT.C

Description:

The Cadence Bed Range from Harvest Healthcare is a high-end profiling bed designed exclusively for the care home market. Its standout feature is the ability to be fully customised to meet the unique needs and preferences of each home. Customers can choose from various options, including width, frame, headboard, footboard, colors, fabrics, side rails, and side skirts, all crafted through in-house upholstery services.

The Cadence Comfort Bed is a stylish, versatile and cost-effective pre-configured option that features a standard head and footboard with a curved design, complete with side skirts. The bed can be used with either continuous side rails or split side rails, which are available as optional extras, providing the flexibility to meet a variety of needs and preferences.

With its elegant design and adaptable features, the Cadence Comfort Bed is a comfortable and stylish choice for any bedroom.

Max Patient Weight:	200kg / 21 Stone
Safe Working Load:	235 kg / 37 Stone
Overall Dimension:	2115 x 1055mm
Adjustable Platform Height:	240mm - 800mm
Mattress Dimensions:	2000 x 900mm (Wide frame available)
Angle of Back Rest:	Up to 70 degrees
Angle of Knee Brake	Up to 45 degrees
Bed Weight:	126 kg

Key Information:

- Trendelenburg
- Auto-Regression
- Extendable
- Standard & Extended Side Rails
- Lockable Handset
- Low Voltage Power Supply
- Emergency Battery Back-Up
- 5 Position Foot Rest
- IPX4 Water Resistance
- Reverse Trendelenburg

Key Information:

- High-end profiling bed
- Fully customisable to meet the specific needs and preferences of each care home
- In-house upholstery services offer a wide range of colours and fabrics for personalisation
- Options to choose width, frame, headboard, footboard, side rails, and side skirts
- Three pre-configured choices – Cadence Comfort, Cadence Select, and Cadence Premium
- 5-year warranty for added peace of mind and long-term support
- Advanced safety and comfort features prioritise resident well-being
- Elevates the care experience with unmatched luxury, style, and functionality

Supplier Contact Details:

t: +44(0)1709 377172

e: enquiries@harvesthealthcare.co.uk

w: www.harvesthealthcare.co.uk

Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products.

Products and data continually undergo development and may therefore deviate from the details given.

Version 1: July 2023