

# Cadence Premium Profiling Bed



**Product Code:**

**CADENCE.PREMIUM**

**Description:**

The Cadence Bed Range from Harvest Healthcare is a high-end profiling bed designed exclusively for the care home market. Its standout feature is the ability to be fully customised to meet the unique needs and preferences of each home. Customers can choose from various options, including width, frame, headboard, footboard, colors, fabrics, side rails, and side skirts, all crafted through in-house upholstery services.

The Cadence Premium bed is a luxurious and elegant pre-configured bed that boasts a curved headboard and low curved footboard, both accented with stylish side skirts. The bed is designed to be used with split side rails, which can be added as optional extras, providing enhanced support and customisation options to meet individual needs.

With its refined design and customisable features, the Cadence Premium bed is a perfect choice for those seeking a high-end and personalised sleeping experience.

<b>Max Patient Weight:</b>	200kg / 21 Stone
<b>Safe Working Load:</b>	235 kg / 37 Stone
<b>Overall Dimension:</b>	2115 x 1055mm
<b>Adjustable Platform Height:</b>	240mm - 800mm
<b>Mattress Dimensions:</b>	2000 x 900mm (Wide frame available)
<b>Angle of Back Rest:</b>	Up to 70 degrees
<b>Angle of Knee Brake</b>	Up to 45 degrees
<b>Bed Weight:</b>	126 kg

**Key Information:**

- Trendelenburg
- Auto-Regression
- Extendable
- Standard & Extended Side Rails
- Lockable Handset
- Low Voltage Power Supply
- Emergency Battery Back-Up
- 5 Position Foot Rest
- IPX4 Water Resistance
- Reverse Trendelenburg

**Key Information:**

- High-end profiling bed
- Fully customisable to meet the specific needs and preferences of each care home
- In-house upholstery services offer a wide range of colours and fabrics for personalisation
- Options to choose width, frame, headboard, footboard, side rails, and side skirts
- Three pre-configured choices – Cadence Comfort, Cadence Select, and Cadence Premium
- 5-year warranty for added peace of mind and long-term support
- Advanced safety and comfort features prioritise resident well-being
- Elevates the care experience with unmatched luxury, style, and functionality

**Supplier Contact Details:**

t: +44(0)1709 377172

e: [enquiries@harvesthealthcare.co.uk](mailto:enquiries@harvesthealthcare.co.uk)

w: [www.harvesthealthcare.co.uk](http://www.harvesthealthcare.co.uk)

Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

**Important Safety Notice**

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products.

Products and data continually undergo development and may therefore deviate from the details given.

**Version 1: July 2023**