

# PivotPro Overbed Table



**Product Code:** HLDWUL120

- Description:**
- Introduces a smart pivot design that securely clamps to the bed frame to address challenges associated with floor beds and bedside tables; traditional bedside tables become impractical when the bed is lowered, and standard overbed tables cannot be used either when the bed is lowered due to clearance under the bed.
  - The PivotPro can be effortlessly pivoted to the side of the bed, transforming into a functional bedside table for lightweight personal items such as remote controls, eyewear, and books.
  - Alternatively, it can be pivoted in front of the person in bed to provide a convenient surface for holding food, drinks, and personal effects.
  - Compatible with other Woburn beds, excluding the Woburn Ultimate.
  - Equipped with removable spacers, it accommodates beds with larger platform tubes, expanding its usage beyond the Woburn range.

**Weight:** 6.0 kg

**Clearance required under bed:** Approx 10mm

**Safe Working Load:** 8.0 kg

**Warranty:**

Subject to the terms of Harvest Healthcare Guarantee and Warranties the PivotPro Overbed Table is guaranteed for one year from date of purchase.

**The Pivot Pro table can be used in 2 positions.**



**Bedside Table**



**Overbed Table**

**Supplier Contact Details:**

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk  
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

**Important Safety Notice**

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

**Version 1: January 2024**