



PARTNERSHIP ACTIVE MATTRESS

**Product Code:**

HR500

Description:

- The Partnership Mattress is a double bed sized mattress with one side foam and one side alternating pressure relieving air cells. The alternating side will provide relief for patients at High Risk of developing pressure ulcers
- Maintains Quality of Life, Enabling Couples to Stay Together
- Visco topped static mattress and foam walls
- Alternating pressure relief
- 10 independent cells
- Foam can be turned due to allow for comfort preference
- Removable cell structure aids infection control
- Low running costs
- Quiet operation
- CPR facility
- Can be used on a divan style bed as well as our 1400mm wide Woburn Ultimate bed
- Complete with quilted vapour permeable water resistant cover and heavy duty base sheet

Mattress Dimensions:

1980 x 1370 x 190mm

Weight:

9kg

Max. Patient Weight:

25 Stone (158kg)

Pressure Relief Risk:

Very High

Pump Technical Specification:

Please see separate pump specification sheet

Cleaning Guidance:

Cover should be wiped down with hot water and mild detergent. If heavily soiled, the cover may be machine washed up to 85°C and tumble dried. For disinfection use sodium hypochlorite solution max 1000ppm. Use no bleach, biological or phenolic-based cleaners. Dry completely before storing.

Warranty

In terms of Harvest Healthcare Guarantee and Warranties the Harvest Partnership mattress is guaranteed for one year from date of purchase.

Exclusions:

1. Removable covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded
3. Components of pumps which are recommended for annual service, are guaranteed for one year.

Supplier Contact Details:

t: +44(0)1709 377172 e: sales@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact.
(Mattress risk categories are to be used as a guide only)

Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 2: Jun 2018