



Lissa Oak

NON-TILTING OVERBED TABLE

LISSA OAK: HMT115LO
OYSTER URBAN OAK: HMT115OU

The versatile Non-Tilting Overbed Table is thoughtfully designed to enhance convenience and safety in any healthcare or home setting. Featuring a rigid tabletop and a sturdy, slim design, it seamlessly integrates into various environments. Its standout one-hand operation and stepless height adjustment (from 760mm to 1000mm) ensure effortless customisation to suit individual user needs and bed requirements, promoting greater independence.

The innovative U-shape chassis allows for smart, space-saving storage directly behind the bedside or bed, optimising room layout and improving accessibility. For added security, a robust metal rail is integrated to prevent personal items, meals, or drinks from falling, providing peace of mind. Built with durable, powder-coated metalwork and designed to be easy to clean, this table supports high hygiene standards and promises long-lasting performance.

HIGHLIGHTS



**One-Hand
Operation**



Easy to Clean



**Stepless Height
Adjustable**



**15 kg High Load
Capacity**

CONTACT US

☎ +44 (0)1709 377172
✉ sales@harvesthealthcare.co.uk

FOLLOW US

in [harvest-healthcare-ltd](#)
x [@HarvestHealthUK](#)

OUR LOCATION

📍 Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

KEY PRODUCT FEATURES

-  **One-Hand Operation:** Designed for ease of use, allowing adjustments with a single hand.
-  **Easy to Clean:** Surfaces are simple to wipe down, supporting high hygiene standards.
-  **Stepless Height Adjustable:** Offers smooth, continuous height adjustment from 760mm to 1000mm.
-  **High Load Capacity:** Capable of supporting a maximum load of up to 15 kg.
-  **Secure Mobility:** 4 individually lockable 75mm castors for stable positioning and easy movement.
-  **Space-Saving Design:** U-shaped chassis allows for efficient storage.
-  **Durable Construction:** Metalwork is powder-coated in RAL 9006 for enhanced durability and a clean finish.

For Further Information:

Scan the QR Code below to visit our website.



TECHNICAL DATA	
Product Code	HMT115LO/HMT115OU
Warranty	1 Year
Weight	15 kg
Dimensions	446 x 900mm
Height Adjustment Range	760-1000mm
Table Top Dimensions	900 x 400 x 190 mm
Height of Wheel Bar from Floor	122.5mm
Castor Diameter	75mm
Safe Working Load (Table Top)	15 kg
Materials	Castor Bracket, Bracket: steel (powder-coated). Table Top: Wood Material (Veneered). Castors, Handles: Plastic.



Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 1: July 2025