

OLYMPIC PRO 185 HOIST



Product Code: HLA002

Description:

- The Olympic Pro 185 is the ideal tool for safe and effortless lifting, repositioning and transferring of residents who are limited in their mobility.

Safe Working Load: 185kg (29 stone)

Chassis Width: Chassis 700mm/1180mm (Closed/Spread)

Length/Height: Total 1290mm/125mm

Width: 700mm

Lifting Range: 1400mm

Maximum Height: 2170mm

Key Features:

- Aesthetic design
- Lightweight aluminium and steel frame
- Soft stop/start
- Lifts from the floor
- Electrically adjustable height, cradle and spreadable carrier frame
- Anti-rotation mechanism
- Cradle with Electrical Positioning System (EPS)
- Lockable double castors
- Emergency stop and emergency lowering
- Hand control with display of battery charge status
- LCD display for service notification, actuator cycles and overloads

Optional Extras:

- In-built electronic scales
- Spare battery
- Wall charging station



Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Olympic Pro Hoist is guaranteed for five years from date of purchase.

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 1: May 2019