

harvest

Together we
transform lives

Cadence

Bed Range



**MEET THE
NEW CADENCE
FLOOR BED**

[harvesthealthcare.co.uk](https://www.harvesthealthcare.co.uk)

Discover

Cadence

The Cadence Bed Range by Harvest


Every detail matters when creating spaces where residents truly feel at home. The Cadence range brings together thoughtful design and genuine comfort – helping you curate a premium care environment with confidence.

One range. Four options. A simple choice.

Now featuring the new Cadence Floor Bed, the collection offers four tiered options to suit every need: Cadence, Cadence Floor, Cadence Upholstered, and Cadence Floor Upholstered.

Match your interior with three classic wooden styles (Comfort, Select, Premium) or three elegant upholstered finishes (Fluted, Regal, Winged).

Discover the Cadence. Discover the difference.







Contents

01 GALLERY

04 CADENCE WOODEN STYLES

07 CADENCE UPHOLSTERED STYLES

10 SPECIFICATIONS

13 OPTIONAL ACCESSORIES

14 WHY CHOOSE HARVEST?

Wooden Cadence



Cadence Comfort



Cadence Select



Cadence Premium

Upholstered Cadence



Cadence Fluted



Cadence Regal



Cadence Winged

The New Cadence Floor



Cadence Floor



Cadence Floor



Cadence Floor





Comfort Style

The Cadence Comfort features high head and footboards with curved corners, with or without channels – giving you the flexibility to tailor the look to your environment.


Compatible with both full-length and split side rails (available as optional extras). Comfort is the only style in the range that supports full-length side rails, making it a versatile and reliable choice for a wide range of needs.

Full Length Siderails only compatible with Comfort Style.

Available In:

-  Lindberg Oak
-  Oyster Urban Oak

Frame: Cadence or Cadence Floor




Headboard: C-Type




Footboard: C-Type



Siderails: 2-Part, 3-Part Split or Full Length



Sideskirts: 2-Part, 3-Part Split or Full Length



Both frames are available in **90cm or 105cm** widths, with a in-built 20cm extension as optional on the Cadence and standard on the Cadence Floor.





Select Style

The Cadence Select pairs a high, square-cornered headboard with a low footboard for a sleek, contemporary profile.


Optional split side rails are available in 2-part or 3-part configurations, each can be paired with optional matching split sideskirts. Without side rails, optional standard sideskirts complete the look – practical and refined without compromise.

Full Length Sideskirts only available without siderails.

Available In:

-  Lindberg Oak
-  Oyster Urban Oak

Frame: Cadence or Cadence Floor



Headboard: V-Type



Footboard: V-Type Low



Siderails: 2-Part or 3-Part Split



Sideskirts: 2-Part, 3-Part Split or Full Length



Both frames are available in **90cm or 105cm** widths, with a in-built 20cm extension as optional on the Cadence and standard on the Cadence Floor.





Premium Style

The Cadence Premium features a curved headboard and low curved footboard for a softer, more refined silhouette.

Optional split siderails in 2-part or 3-part configurations can be added for enhanced support, each paired with optional matching split sideskirts. Without siderails, optional full length sideskirts can complete the look – elegant and tailored to suit.

Full Length Sideskirts only available without siderails.

Available In:

-  Lindberg Oak
-  Oyster Urban Oak

Frame: Cadence or Cadence Floor



Headboard: Curved




Footboard: Curved Low



Siderails: 2-Part or 3-Part Split



Sideskirts: 2-Part, 3-Part Split or Full Length



Both frames are available in **90cm or 105cm** widths, with a in-built 20cm extension as optional on the Cadence and standard on the Cadence Floor.






7
YEARS
WARRANTY*

Fluted Style

The Cadence Fluted Style features a high-angled headboard with refined vertical channels and a low, clean footboard for a sophisticated, contemporary silhouette.

Choose from a curated fabric selection to seamlessly match your room décor. Optional split side rails in 2-part or 3-part configurations can be added, each paired with optional matching split sideskirts. Without side rails, full length sideskirts can complete the look – architectural and elegant.

Available In:

-  Compass Denim
-  Compass Canvas
-  Compass Silver

Frame: Cadence or Cadence Floor



Headboard: Fluted





Footboard: Upholstered Low



Siderails: 2-Part or 3-Part Split



Sideskirts: 2-Part, 3-Part Split or Full Length



Both frames are available in **90cm or 105cm** widths, with a in-built 20cm extension as optional on the Cadence and standard on the Cadence Floor.






7
YEARS
WARRANTY*

Regal Style

The Cadence Regal Style features an elegant arched headboard with a refined crafted accent and a low footboard, bringing timeless sophistication and graceful curves to the care environment.

Choose from a range of fabric colours to complement your décor. Optional split siderails in 2-part or 3-part configurations can be added, each paired with optional matching split sideskirts. Without siderails, full length sideskirts complete the look – classic and distinguished.

Available In:

-  Compass Denim
-  Compass Canvas
-  Compass Silver

Frame: Cadence or Cadence Floor



Headboard: Regal





Footboard: Upholstered Low



Siderails: 2-Part or 3-Part Split



Sideskirts: 2-Part, 3-Part Split or Full Length



Both frames are available in **90cm or 105cm** widths, with a in-built 20cm extension as optional on the Cadence and standard on the Cadence Floor.






7
YEARS
WARRANTY*

Winged Style


The Cadence Winged features statement winged sides on the headboard and a low footboard, adding depth and a sense of warmth to create a cosy, enveloping feel.

Choose from a variety of fabric colours to match your décor. Optional splitside rails in 2-part or 3-part configurations can be added, each paired with optional matching split sideskirts. Without siderails, full length sideskirts complete the look – distinctive and inviting.

Available In:

-  Compass Denim
-  Compass Canvas
-  Compass Silver

Frame: Cadence or Cadence Floor



Headboard: Winged




Footboard: Upholstered Low



Siderails: 2-Part or 3-Part Split



Sideskirts: 2-Part, 3-Part Split or Full Length



Both frames are available in **90cm or 105cm** widths, with a in-built 20cm extension as optional on the Cadence and standard on the Cadence Floor.











Cadence Specs

The Cadence brings together comfort, ergonomics, and considered design in one complete care bed range.

The following pages outline key specifications across the range – from dimensions and weight capacity to mobility features and positioning capabilities including Trendelenburg, auto-regression, and fully adjustable height and angle settings.

Key Features

 Warranty	*7 Years
 Trendelenburg	Yes
 Auto-Regression	Yes
 Max User Weight	200kg 31.5st
 Safe Working Load	235kg 37st
 Dimensions (900)	2155 x 1055mm 7 x 3.5ft
 Height Range	240 - 800mm 0.8 x 2.6ft
 Bed Weight	126kg 19.8st

FOR MORE INFORMATION, VISIT OUR WEBSITE

*7-Year Warranty on Frame. Plus, 5-year warranty on woodwork and electrics.

7
YEARS
WARRANTY*



Cadence Specs Cont.

Combining refined aesthetics with a comprehensive range of features to meet evolving resident needs while promoting comfort and independence.



Backrest
0 - 70°



Thighrest
0 - 45°



Trendelenburg
Included



Height Range
240 - 800mm

Additional Features

- ✓ Optional Platform Extension
- ✓ Standard & Extended Siderails
- ✓ 70 Degree Back Rest Angle
- ✓ 45 Degree Thigh Rest Angle
- ✓ Lockable Handset
- ✓ Low Voltage Power Supply
- ✓ Emergency Battery Back-Up
- ✓ Manual & Electric Calf Adjustment
- ✓ IPX4 Ingress Protection
- ✓ Reverse Trendelenburg
- ✓ Optional Central Brake
- ✓ Single Function Chair Position
- ✓ Integral Backrest Autoregression

FOR MORE INFORMATION, VISIT OUR WEBSITE

*7-Year Warranty on Frame. Plus, 5-year warranty on
woodwork and electrics.



Cadence Floor Specs

The following page outlines the key specifications of the Cadence Floor – from dimensions and platform height range to weight capacity, mobility features, and advanced positioning capabilities.



Backrest
0 - 80°



Thighrest
0 - 35°



Trendelenburg
Included



Height Range
110 - 800mm

Warranty	*7 Years
Trendelenburg	Yes
Auto-Regression	Yes
Max User Weight	170kg 26.8st
Safe Working Load	205kg 32.3st
Dimensions	2100 x 1050mm 6.9 x 3.4ft
Height Range	110 - 800mm 0.4 x 2.6ft
Bed Weight	142kg 22.3st

- ✓ Two-Stage Lowering
- ✓ Integrated Platform Extension
- ✓ Standard & Extended Siderails
- ✓ 80 Degree Back Rest Angle
- ✓ 35 Degree Thigh Rest Angle
- ✓ + All Other Additional Features
as per the standard Cadence bed.

FOR MORE INFORMATION, VISIT OUR WEBSITE

*7-Year Warranty on Frame. Plus, 5-year warranty on
woodwork and electrics.



Optional Accessories

Every environment is different, and every resident's needs are unique. Our range of optional accessories lets you tailor the Cadence to suit – from practical additions like the integral bed extension, lifting pole, and grab handle, to comfort-focused extras like the siderail cushion. Under-bed and over-bed lighting options provide gentle, low-level illumination for safer movement during the night.

Bed Extension



Under Bed Light



Over Bed Light




Lifting Pole



Grab Handle



Siderail Bumper



FOR MORE INFORMATION, VISIT OUR WEBSITE

Why Choose Harvest?

We believe great care environments start with great partnerships. From your first conversation to years down the line, we're with you.

Honest, friendly service and a responsive team that works around your schedule – including delivery, installation, and setup on a date that suits you. Our nationwide engineering team provides lifetime servicing, maintenance backup, and readily available parts to keep everything running smoothly. We also offer MO training for bed care and troubleshooting, so your team always feels confident.

Our in-house upholstery and design team means every detail is considered and quality controlled from start to finish – giving you unrivalled design flexibility and a true one-stop-shop experience we're proud to stand behind.

Every Cadence bed is built to exceptional standards and backed by a 7-year warranty*, with full asset register provision as standard.





harvest
Together we
transform lives

GET IN TOUCH



+44 (0)1709 377172



[harvesthealthcare.co.uk](https://www.harvesthealthcare.co.uk)



sales@harvesthealthcare.co.uk



Sheaf House
Bradmarsh Business Park
Bradmarsh Way
Rotherham
S60 1BW



[harvest-healthcare-ltd](https://www.linkedin.com/company/harvest-healthcare-ltd)



[@HarvestHealthUK](https://twitter.com/HarvestHealthUK)



[Harvest Healthcare](https://www.youtube.com/channel/UC...)



Visit the
Cadence
Webpage