



Cadence Profiling Bed

The Cadence Bed Range reimagines resident comfort and elevates your care environment. Each bed blends advanced technology with premium comfort, creating a space that resonates with residents and families and contributes to your home's overall value proposition. A standout feature is the ability to be fully customised to meet the unique needs and preferences of each home.

With a 7-year warranty on the bed frame and advanced safety and comfort features, this bed range ensures a great level of care, comfort, and style for residents, promoting a truly tailored care experience. The Cadence is ideal for long-term care environments where resident comfort, safety, and individualised care are paramount. Its versatility, customisation options, and focus on resident well-being make it a valuable asset in any care setting.

Key Product Features



Premium, Customisable Design: Creates an elevated, homely, and sophisticated care environment, moving beyond a clinical appearance.



Customisable Upholstery: Personalised beds with a variety of elegant headboard and footboard options to improve the look and feel of your rooms.



Three Style Choices: Cadence Comfort, Cadence Select, and Cadence Premium to suit a range of different environments and requirements.

HIGHLIGHTS



***7-Year Warranty on Frame.**



Unmatched Customisation



Long-Term Value & Durability



Designed to Elevate Your Spaces

CONTACT US

+44 (0)1709 377172
sales@harvesthealthcare.co.uk

FOLLOW US

harvest-healthcare-ltd
@HarvestHealthUK

OUR LOCATION

Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

Additional Features

- Standard Trendelenburg and Reverse Trendelenburg.
- Standard Auto Regression: improved positioning and reduced shear/friction.
- Optional Inbuilt Extendable Platform
- Optional Standard & Extended Side Rails
- Low Voltage Power Supply
- Lockable Handset (all functions or Trendelenburg only).
- Emergency Battery Back-Up

For further information:

Scan the QR Code below to visit our website.



TECHNICAL DATA

| | |
|------------------------------------|--|
| Product Code | For product codes and ordering, get in touch with our team. |
| Warranty | 7-Year Warranty on Frame. Plus, 5-year warranty on woodwork and electrics. |
| Total Bed Weight | 126 kg |
| Platform Dimensions | 2000 x 900mm |
| Adjustable Height | 240mm-800mm |
| Safe Working Load | 235 kg / 37 Stone |
| Maximum User Weight | 200kg / 31 Stone |
| Overall Dimensions | 2155 x 1055mm |
| Angle of Back Rest | Up to 70° |
| Angle of Knee Rest | Up to 45° |
| Motor Unit Protection Class | IPX4 |
| Medical Product | Class 1 |

ACCESSORIES

| | |
|----------------------------------|---------|
| Universal Safe Roll | HFB003 |
| Grab Handle | HLB685 |
| Grab Handle Safety Sleeve | HLBGHS |
| Lifting Pole | HLDW582 |

Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 3: Feb 2026