

harvest

Together we
transform lives



**ProForm
Range**

Contents

	Page
Introduction	3
Prism Healthcare	4
Multi C-Air	6
Madison	8
ProForm Porter	10
ProForm Position	12
ProForm Contour	14
Boston	16
Cushion Guide	18
Fabrics	20
Warranty	21



Introduction

You already know us for pressure care, profiling beds, and moving and handling equipment. Now we can help you with seating too.

Repose Furniture have been designing and manufacturing bespoke healthcare, hospital, and riser recliner chairs since 1999. Like Harvest, they're part of the Prism Healthcare Group, and they share the same approach we do, building each chair around the individual, with the person's needs and clinical requirements at the centre of the design.

For you, that means one more part of your care environment covered by a team you already trust. Same standards, same commitment to quality and value, same people to call.


Together, we transform lives.



Prism Healthcare

Prism Healthcare is a leading integrated solutions provider dedicated to transforming lives in the paediatric, home, and long-term care sectors. As a collection of companies, we offer a comprehensive range of products and services designed to enhance care quality, improve well-being, and ultimately transform lives.

When health or mobility declines, the options to live independently can become limited. We take pride in offering a range of products and services tailored to empower individuals with limited mobility, helping them regain their independence and enhance their quality of life.



As a UK manufacturer, Prism Healthcare offers affordable and clinically effective products. Our portfolio includes equipment for safe patient handling, bathing, accessibility, speciality seating and pressure area care. These solutions empower healthcare professionals and enable individuals to live life to the fullest with dignity, whether residing in their own homes or long-term care facilities.

Prism Healthcare is not only a leading manufacturer and supplier of a vast portfolio of specialist equipment, but we also go above and beyond to support the care community. We have a rich history and understand that providing exceptional care requires more than just products.

Together, we transform lives.





Multi C-Air

The success of the Multi C-Air lies in its complete adaptability. This riser recliner puts the healthcare team, in the home or the ward, in ultimate control, providing a full range of options to tailor the chair precisely for each individual’s comfort and care. With a range of pressure and postural management cushions, this model is ideal for deteriorating conditions.

Key Features:

- Choice of interchangeable postural management back cushions
- Choice of interchangeable pressure management seat cushions
- Choice of rise and recline mechanisms with the option of VAL technology
- Seat and leg rest generously padded with memory foam
- Heavy duty glides built-in as standard
- Simple to use handset
- Optional seat depth adjustment
- Optional seat height adjustment
- Optional full length arm covers
- Optional safety features
- Range of additional headrests.

Size Guide

The Multi C-Air has a weight limit of 127kg (20 stone) as standard on the small model and 159kg (25 stone) on the medium and large.

	Small	Medium	Large
Seat Height	16" (406mm)	18" (457mm)	20" (508mm)
Seat Width	18" (457mm)	20" (508mm)	22" (559mm)
Seat Depth	18" (457mm)	20" (508mm)	20" (508mm)
Back Height	28/29" (711/737mm)	28/29" (711/737mm)	28/29" (711/737mm)
Arm Height	7" (178mm)	7" (178mm)	7" (178mm)





Madison

The Madison has all the benefits of a healthcare riser recliner chair, with the added versatility of the large lockable castors for moving clients around. The Dual Tilt-in-Space mechanism maximises the pressure and posture management options for the user with independent movements, as well as a rise function for ease of standing.

Key Features:

- Dual Tilt-in-Space mechanism with independent movements
- Choice of interchangeable postural management back cushions
- Choice of interchangeable pressure management seat cushions
- A single chaise cover to reduce pinch points and prevent shear
- Heavy-duty lockable castors
- Handset storage loop on the arm
- Heavy-duty push handle with rubber grips for portering
- Optional seat depth adjustment
- Optional sliding footplate
- Optional full length arm covers and head cover
- Optional safety features and headrests
- Available in petite sizes and as an express chair.

Size Guide

The Madison has a weight limit of 127kg (20 stone) as standard.

	Small	Medium	Large
Seat Height	19" (483mm)	19" (483mm)	19" (483mm)
Seat Width	18" (457mm)	20" (508mm)	22" (559mm)
Seat Depth	20" (508mm)	20" (508mm)	20" (508mm)
Back Height	28" (711mm)	28" (711mm)	28" (711mm)
Arm Height	6" (152mm)	6" (152mm)	6" (152mm)





ProForm Porter

Every aspect of the ProForm Porter healthcare chair has been researched and engineered to provide a feeling of comfort, safety, and security. From the high arm design and backrest recline to the seat depth and angle adjustments, each element of the design is tailored to support conditions associated with involuntary or ataxic movement.

Key Features:

- Manual mechanism with independent backrest recline
- Choice of interchangeable postural management back cushions
- Choice of interchangeable pressure management seat cushions
- Seat angle adjustment
- Seat depth adjustment
- A range of healthcare fabric options
- Sliding footplate
- Heavy-duty push handle with rubber grips
- Heavy-duty lockable castors
- D-rings
- Optional leg support block
- Optional connecting straps for a leg support block.

Size Guide

The ProForm Porter has a weight limit of 127kg (20 stone) as standard.

	Small	Medium	Large
Seat Height	16" (406mm)	18" (457mm)	21" (533mm)
Seat Width	17" (432mm)	19" (483mm)	21" (533mm)
Seat Depth	17-21" (432-533mm)	17-21" (432-533mm)	17-21" (432-533mm)
Back Height	28" (711mm)	28" (711mm)	28" (711mm)
Arm Height	6.5" (165mm)	6.5" (165mm)	6.5" (165mm)





ProForm Position

The ProForm Position is a simple, but multi-adjustable seating solution. Changing user requirements and the need for long-term sustainability across a number of environments have been a key part of the ProForm Position’s design. The ability to adjust, replace, and reuse means this model delivers excellent long term user and clinical benefits.

Key Features:

- Available in either an electric or manual model
- Independent backrest recline, leg rest elevation and Tilt-in-Space
- Choice of interchangeable postural management back cushions
- Choice of interchangeable pressure management seat cushions
- Heavy-duty lockable castors
- Heavy-duty push handle with rubber grips
- Tool-free seat depth adjustment
- Tool-free height and angle adjustable footplate
- Seat width adjustment
- Back height adjustment
- An optional foam footpad, available in 3 sizes, and optional directional castor.

Size Guide

The ProForm Position has a weight limit of 160kg (25 stone) as standard.

Seat Height	24.5" (622mm)
Seat Width	16.5-20" (420-508mm)
Seat Depth	15-21" (381-533mm)
Back Height	27-29" (685-736mm)
Arm Height	6.5" (165mm)





ProForm Contour

The ProForm Contour is a multi-adjustable seating solution that has been designed to help users at risk of pressure injuries and falls. We are conscious of changing user requirements and the need for long term sustainability. The ability to adjust, replace and reuse, along with clinical features, has been at the core of the ProForm Contour design.

Key Features:

- Available in either an electric or manual model
- Independent backrest recline, leg rest elevation and Tilt-in-Space
- Choice of interchangeable postural management back cushions
- Choice of interchangeable pressure management seat cushions
- Heavy-duty lockable castors
- Ergonomic back and wing handles
- Tool-free seat depth adjustment
- Tool-free height and angle adjustable footplate
- Seat width adjustment
- Back and arm height adjustment
- Height, depth, and angle adjustable headrest.

Size Guide

The ProForm Contour has a weight limit of 160kg (25 stone) as standard.

Seat Height	24.5/26" (622/660mm)
Seat Width	17.5-22" (445-558mm)
Seat Depth	15-21" (381-533mm)
Back Height	29.5-33.5" (750-850mm)
Arm Height	6-10.5" (152-266mm)





Boston

The Boston is an extremely comfortable healthcare chair that can be used as a general porter chair, but is equally suited to conditions requiring a 24-hour care programme. The independent movement of the Tilt-in-Space, backrest recline and leg rest elevation functions enables a whole range of postural positions for high level pressure management care.

Key Features:

- Available in either an electric or manual model
- Independent backrest recline, leg rest elevation and Tilt-in-Space
- Choice of interchangeable postural management back cushions
- Choice of interchangeable pressure management seat cushions
- Heavy-duty lockable castors
- Seat and leg rest generously padded with memory foam
- D-rings
- Optional seat depth adjustment
- Optional seat width adjustment
- Optional height and angle adjustable footplate
- A range of other optional accessories.

Size Guide


The Boston has a weight limit of 159kg (25 stone) as standard.

	Small	Medium	Large
Seat Height	25" (635mm)	25" (635mm)	25" (635mm)
Seat Width	18" (457mm)	20" (508mm)	22" (559mm)
Seat Depth	20" (508mm)	20" (508mm)	20" (508mm)
Back Height	30" (762mm)	30" (762mm)	30" (762mm)
Arm Height	7" (178mm)	7" (178mm)	7" (178mm)



Cushion Guide

Back Cushions




Our range of removable and interchangeable back cushions are designed to promote optional posture and assist with alignment. These cushions aid in reducing pressure and providing clinical support for a wide range of needs.

Each cushion is quick and simple to change, so the seating solution can adapt to changing needs. For more information on our available cushions and model compatibility, please visit our website.


Seat Cushions

Viscoform




A high-density foam base that offers a supportive and firm base, with a top layer of memory foam to allow immersion into the cushion, providing stability and a large surface area to reduce pressure levels.

Liquiform



A high density foam that offers a supportive and firm seat base. It's topped with fluid liquid gel in a soft silicone cover which provides high pressure relief and helps dissipate heat, while adapting to movement.

Dynaform



A dynamic alternating air cushion which constantly changes the air levels in its air cells to recreate normal body movements. This provides constant pressure redistribution to prevent pressure build up.

The removable and interchangeable seat cushions we offer have been developed to provide the utmost comfort with varying levels of pressure management. Whether for everyday seating, high-risk pressure care or pressure redistribution, there is a cushion solution tailored to individual's requirements.

Each cushion is quick and simple to change, so the seating solution can adapt to changing needs. For more information on our available cushions and model compatibility, please visit our website.



Fabrics

All our healthcare fabrics are Waterproof, Anti-MRSA, Anti-bacterial, Anti-fungal, Crib 5 and bleach cleanable. Our cushions and cushion covers are also all removable for ease of cleaning.

For our full range of healthcare fabrics, please visit our website.



Warranty

Repose offer a best-in-class 5 year warranty as standard on all of their products. This offers peace of mind on the structural frame, mechanical and electrical components, and foam and cushioning.

In association with Castelan group, they also offer a 5 year Care Plan, to protect against a comprehensive range of household accidents that result in damage or staining of the product. With the Care Plan, you are covered against:



Structural Defects

- ✓ Breakage or separation of frame components
- ✓ Breaking or bending of metal mechanisms or other metal components
- ✓ Defective mechanical and electrical recliner mechanisms
- ✓ Broken zips and castors
- ✓ Peeling of leather and peeling/lifting of veneer
- ✓ Electronic equipment, such as USB's
- ✓ Cushion interiors and webbing

The Benefits

- ✓ 5 year peace of mind
- ✓ 24/7 online claim service
- ✓ No excess charges
- ✓ National network of professional technicians
- ✓ No call out fees

Stains

- ✓ Tea, coffee, hot chocolate and milk
- ✓ Shoe polish
- ✓ Wine, beer and spirits
- ✓ Cosmetics
- ✓ Tomato ketchup
- ✓ Soap products
- ✓ Paint
- ✓ Blackcurrant and orange juice
- ✓ Tar
- ✓ Human and animal bodily fluids
- ✓ Cola and other fizzy drinks
- ✓ Mineral oil and glue
- ✓ Curry and pizza
- ✓ Grease from foodstuffs
- ✓ Dye transfer
- ✓ Chocolate
- ✓ Corrosive substances
- ✓ Bleach

Accidental Damage

- ✓ Rips
- ✓ Tears
- ✓ Burns
- ✓ Scuffs
- ✓ Scratches
- ✓ Punctures
- ✓ Dents and chips
- ✓ Pet damage
- ✓ Transit damage
- ✓ Breakage of frame components





harvest
Together we
transform lives

Harvest Healthcare
Sheaf House
Bradmarsh Business Park
Bradmarsh Way
Rotherham
S60 1BW

Tel: +44 (0)1709 377172
Email: sales@harvesthealthcare.co.uk



www.harvesthealthcare.co.uk



**ProForm
Range**

Registered Company Number: 03727157 © Copyright 2026 Repose Furniture Ltd reserves the right to modify or change the specification without prior notification. For terms and conditions of sale and or agreement please refer to your supplier. This brochure is produced as a guide only and colours are accurate only as far as printing tolerances allow.

REP-PRO-RH-0726