



Balmoral II Active Mattress

The Balmoral II is Harvest's premium pressure-reducing active mattress, designed for high-risk users and effective in managing up to Category 4 pressure ulcers. It features 27 PU cells (23 alternating, 4 static) on a 50mm foam base, ensuring optimal pressure redistribution and comfort. It is compatible with Harvest 250 and Harvest 3 Pumps.

Key Product Features

-  **27 PU Alternating Cells:** 23 active and 4 static cells providing effective pressure relief for very high-risk users.
-  **50mm Foam Sub-Mattress:** Enhanced comfort, stability and safety in the event of a power outage.
-  **Versatile Pump Modes:** Alternating, Constant Low Pressure, Max Inflate and Seated modes for flexible care.
-  **Easy-Clean Protective Cover:** 4-way stretch, vapour-permeable, antibacterial PU cover with Crib 5 fire rating.

HIGHLIGHTS

-  **Dual Therapy**
1 in 2 Alt. & CLP
-  **Variable Comfort**
Settings
-  **Max Inflate**
Function
-  **Seated Mode**
-  **No Minimum**
Weight Limit

CONTACT US

+44 (0)1709 377172
sales@harvesthealthcare.co.uk

FOLLOW US

harvest-healthcare-ltd
@HarvestHealthUK

OUR LOCATION

Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

Additional Features

- Pump auto lockout prevents unintentional use.
- Dial CPR reduces risk of unintentional activation.
- Pump features audible and visual alarms.
- Transport cap maintains air pressure for approx. 24 hours.
- All cells are interchangeable for ease of repair and cleaning.
- Integrated pillow section for patient stability.
- Lower profile mattress increases bedrail compliance.
- Compatible with Harvest 250 and Harvest 3 pumps.

For further information:

Scan the QR Code below to visit our website.



TECHNICAL DATA

Product Code	HAMBALM003
Warranty	3 Years
Mattress Dimensions (mm)	2000(L) x 900(W) x 150(H)
Mattress Weight	10 kg
Max. User Weight (250 Pump)	200 kg / 31 st
Max. User Weight (3 Pump)	254 kg / 40 st
Min. User Weight	N/A
Pressure Risk	Very High
Alternating Cycle	1 in 2 over 10 mins
Number of Cells	27 (23 alternating, 4 static)
Functions	Alternating, CLP, Max Inflate, Seated
Cover	4-way stretch PU, Crib 5 rated
Medical Product	Class 1

PRODUCT CODES

HAMBALM003.P250	HAMBALM003.P250
HAMBALM003.P3	HAMBALM003.P3
HAMBALM003.HPU250	HAMBALM003.HPU250
HAMBALM003.HPU3	HAMBALM003.HPU3

Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 3: Feb 2026