



Blenheim Active Cushion

The Blenheim Active Cushion and 250 Pump provides a low-profile, alternating and low-pressure cushion system of 9 full air cells. Designed to support the prevention and management of pressure ulcers for patients at very high risk, including existing tissue damage up to and including Category 4. The low-profile design improves seating posture and stability and aids transfers, while achieving optimum pressure relief and comfort when seated.

KEY PRODUCT FEATURES

- ✓ **9 Full Air Cells:** Alternating therapy in a 1 in 2 over a 10 minute cycle for effective pressure care.
- ✓ **Low Profile Design:** Improves seating posture and stability while aiding patient transfers.
- ✓ **Versatile Pump Modes:** Alternating, Constant Low Pressure, Max Inflate and Seated modes.
- ✓ **Anti-Slip Base:** Maintains patient positioning, with a 4-way stretch welded PU cover.

HIGHLIGHTS

- ★ **9 Cell Alternating Therapy**
- ★ **Very High Risk Pressure Relief**
- ★ **152kg Max User Weight**
- ★ **Pump Auto Lock Out**

CONTACT US

☎ +44 (0)1709 377172
✉ sales@harvesthealthcare.co.uk

FOLLOW US

in harvest-healthcare-ltd
x @HarvestHealthUK

OUR LOCATION

📍 Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

ADDITIONAL FEATURES

- Max Inflate with auto return to alternating for patient safety.
- Seated mode to recognise when the cushion is connected.
- Pump maintains set pressures to reduce carer intervention.
- No minimum user weight limit.
- 4-way stretch, vapour permeable, anti-bacterial welded PU cover with 360° covered zip.
- Pump Auto Lock Out function to reduce risk of unintentional use.
- Fully modular system for ease of repair and cleaning.

For further information:

Scan the QR Code below to visit our website.



TECHNICAL DATA	
Product Code	HAMBALM004.HPU250
Warranty	3 Years
Cushion Dimensions	2000 x 900 x 150 mm
Pump Dimensions	275 x 116 x 175 mm
Max. User Weight	200 kg / 31 st
Min. User Weight	N/A
Total System Weight	11.7 kg
Number of Cells	27 (23 alternating, 4 static)
Pressure Risk	50mm air
Cycle Time	Very High
Medical Product	Class 1

PRODUCT CODES	
HACBLEN001.HP250	Blenheim Active Cushion & 250 Pump

Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 3: June 2026