



Community Static Mattress

The Community Static Mattress is a cost-effective, comfort mattress for use on profiling beds. Built with a high-resilience, fire-retardant foam that allows for better pressure distribution, it is very comfortable for the user and suitable for patients at medium risk of developing a pressure ulcer. Shaping enables profiling of the mattress, and the ability to turn the mattress extends its life. Additional sizes and vandal proof versions are also available.

KEY PRODUCT FEATURES

- ✓ **High Resilience Foam Base:** Fire retardant foam giving better pressure distribution and comfort.
- ✓ **'V' Cuts for Profiling:** Shaping enables the mattress to profile with the bed platform.
- ✓ **Turnable for Longer Life:** The ability to turn the mattress extends its usable service life.
- ✓ **Easy-Clean PU Cover:** 4-way stretch, vapour permeable, water-resistant PU cover.

HIGHLIGHTS

- ★ **Medium Risk Pressure Relief**
- ★ **25 Stone Max Patient Weight**
- ★ **'V' Cuts to Assist Profiling**
- ★ **Cost Effective Quality Build**

CONTACT US

☎ +44 (0)1709 377172
✉ sales@harvesthealthcare.co.uk

FOLLOW US

in harvest-healthcare-ltd
x @HarvestHealthUK

OUR LOCATION

📍 Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

ADDITIONAL FEATURES

- Suitable for patients at medium risk of developing a pressure ulcer.
- Cost-effective solution without compromising on quality.
- Designed for use on electric profiling bed frames.
- Available in a castellated version for increased pressure reduction.
- Additional sizes available, including 1200mm wide (HSMCOM.12).
- Vandal-proof versions available for challenging environments.
- Cover can be wiped down or machine-washed; helps maintain infection control.

For further information:

Scan the QR Code below to visit our website.



TECHNICAL DATA

Product Code	HSMCOM
Mattress Dimensions	2000 x 900 x 125 mm
Max. Patient Weight	159 kg / 25 st
Pressure Risk	Medium
Foam Construction	High resilience, fire retardant foam
Cover	4-way stretch PU, vapour permeable
Profiling	'V' cuts to assist profiling
Turning Required	Turnable to extend mattress life
Medical Product	Class 1

PRODUCT CODES

HSMCOM	Community Static Mattress – 2000 x 900 x 125 mm
HSMCOM.12	Community Static Mattress – 1200mm Wide
HSMCOM.CAS	Community Static Mattress – Castellated

Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 4: July 2026