



Duke Extra 2.0 Active Mattress

The Duke Extra 2.0 is an alternating pressure-relieving mattress system for the prevention and treatment of pressure ulcers. Recommended for patients at very high risk, it treats up to stage 4 pressure ulcers. Made up of 25 alternating cells and 4 static pillow cells, with a 50 stone maximum patient weight, it suits larger patients. Edge-to-edge alternating cells give pressure reduction across the full width, and it fits wider beds like the Community 1200 and Ultimate 1200.

KEY PRODUCT FEATURES

- ✓ **29 PU Cells:** 25 active and 4 static pillow cells with edge to edge alternating pressure care.
- ✓ **50 Stone Capacity:** Designed for larger patients, fitting wider Community 1200 and Ultimate 1200 beds.
- ✓ **50mm Air Sub Mattress:** Provides a stable surface and reduces the risk of bottoming out.
- ✓ **Selectable Cycle Times:** 1 in 2 as default with 10, 15, 20 and 25 minute cycle options.

HIGHLIGHTS

- ★ **50 Stone Max Patient Weight**
- ★ **50mm Air Sub Mattress**
- ★ **Very High Risk up to Stage 4**
- ★ **Fits 1200mm Wide Beds**

CONTACT US

☎ +44 (0)1709 377172
✉ sales@harvesthealthcare.co.uk

FOLLOW US

in harvest-healthcare-ltd
x @HarvestHealthUK

OUR LOCATION

📍 Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

ADDITIONAL FEATURES

- Static function for when a firm, stable surface is needed.
- Max Inflate with auto return to alternating for patient safety.
- Seated mode to support the patient when the backrest is elevated.
- Adjustable pressures to manage individual comfort needs.
- Pump Auto Lock Out function to reduce risk of unintentional use.
- 4-way stretch, vapour-permeable, anti-bacterial PU cover with 360° covered zip.
- Foam base to provide a stable surface and reduce bottoming out.
- Fully modular system for ease of repair and cleaning.

For further information:

Scan the QR Code below to visit our website.



TECHNICAL DATA

Product Code	HAMDUKE.12.HPU3
Warranty	3 Years
Mattress Dimensions	2000 x 1200 x 200 mm
Max. Patient Weight	317 kg / 50 st
Min. User Weight	N/A
Number of Cells	29(25 active, 4 static)
Sub Mattress	50mm air / foam base
Pressure Risk	Very High (up to Stage 4)
Cycle Options	10 / 15 / 20 / 25 mins
Pump	Harvest 3
Medical Product	Class 1

PRODUCT CODES

HAMDUKE.12.HPU3	Duke Extra 2.0 Mattress with Harvest 3 Pump
-----------------	---

Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 3: June 2026