



Hampton II Active Mattress

The Hampton II Active Mattress and 250 Pump is a full mattress replacement system designed to prevent and manage pressure ulcers for patients at very high risk, including those with existing tissue damage up to and including Category 4. It consists of 20 individual cells, with a cell-on-cell construction that provides a stable support base even in the event of a power failure.

The compact digital pump offers automatic pressure regulation and a range of functions, including:

Alternating, Constant Low Pressure, Max Inflate, and Seated modes. This high-specification system is durable, reliable, and easy to use, making it ideal for Long-Term Care and Home Care environments.

The Hampton II Mattress is available in 900mm and 1050mm widths and 2000mm and 2200mm lengths and compatible with all Harvest pumps.

HIGHLIGHTS



**Very High Risk
Pressure Relief**



**No Minimum User
Weight**



**31 Stone
Max User Weight**



**3-Year
Warranty**

CONTACT US

+44 (0)1709 377172
sales@harvesthealthcare.co.uk

FOLLOW US

harvest-healthcare-ltd
@HarvestHealthUK

OUR LOCATION

Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

KEY PRODUCT FEATURES

- ✓ **Cell on Cell Construction:** The lower cell remains inflated to provide a supportive base, ensuring stability and safety in the event of a power failure.
- ✓ **Nylon with TPU Lining Cells:** These cells are designed to reduce the "rustle" noise often associated with active mattresses, enhancing resident comfort and sleep quality.
- ✓ **Versatile Pump Functions:** Provides alternating therapy, Constant Low Pressure, Max Inflate, and Seated modes to meet diverse patient needs.
- ✓ **User-Centric Design:** Features a stable pillow section to enhance patient comfort, stability, and sleep quality.
- ✓ **PU Cover:** A breathable, water-resistant, multi-stretch, and antibacterial cover that is easy to clean and maintain, supporting high hygiene standards.

For further information:

Scan the QR Code below to visit our website.



TECHNICAL DATA	
Product Code	HAMHAMP005.HPU250 (+ variants)
Description	Hampton II Mattress & Harvest 250 Pump
Mattress Dimensions - HAMPHAMP005	2000(L)mm x 900(W)mm x 150(H)mm – Standard Version
HAMPHAMP003	2200(L)mm x 900(W)mm x 200(H)mm – Longer Version
HAMPHAMP1050	2000(L)mm x 1050(W)mm x 200(H)mm – Wider Version
Harvest 250 Pump Dimensions	275(L) x 116(W) x 175(H)mm
Max. User Weight	Max. User Weight 200kg/31st
Min. User Weight	No Minimum User Weight
Mattress Weight	8kg
Total System Weight	10kg
Pressure Risk	Very High Risk
Medical Product	Class II

PRODUCT CODES	
HAMHAMP005.HPU1/HAMHAMP1050.HPU1	Hampton II Mattress System with Harvest 1 Pump (25 Stone / 158kg)
HAMHAMP005.HPU250/HAMHAMP003.HPU250/ HAMHAMP1050.HPU250	Hampton II Mattress with Harvest 250 Pum (31 Stone / 200kg)
HAMHAMP005.HPU3/HAMHAMP1050.HPU3	Hampton II Mattress System with Harvest 3 Pump (40 Stone / 254kg)

Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 4: July 2026